
APEX Series APEX-Serie

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English – page 3

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You have purchased a quality precision weighing instrument that requires handling with care.

Please read these operating instructions carefully before operating your new Denver Instrument balance.

 Calibrate your balance using a reference weight of the appropriate accuracy class. A measuring instrument is only as accurate as the standard that was used as a reference in calibration. For assistance in the selection of reference weights, please contact the factory.

 **Caution!**

The operator shall be responsible for any modifications to Denver Instrument equipment and for any connections of cables or equipment not supplied by the Denver Instrument GmbH and must check and, if necessary, correct these modifications and connections.

 **Warning**

Never lift balance by the draft shield, as this may cause damage to internal mechanisms.

Always lift and transport the balance by its base, including when removing packing materials.

Contents

Specifications	5	Counting	13
Intended Use	6	Weighing in Percent	15
Warning and Safety Information	6	Configuring the Balance	16
Getting Started	7	Device Parameters (Overview)	18
Warranty	7	Troubleshooting	20
Unpacking the Balance	7	RS-232C Interface	21
Installation Instructions	7	Data Input Format	21
Preparing the		Hardware Handshake	22
Weighing Chamber	8	Pin Assignment Chart	22
Connecting the Balance		Care and Maintenance	23
to AC Power	8	Safety Inspection	23
Leveling the Balance	8	Accessories	24
Warmup Time	8	Declaration of Conformity	25
Operating the Balance	11		
Taring the Balance	11		
Toggling between Weight Units	11		
Selecting a Weight Unit	11		
Calibration/Adjustment	12		

Specifications

Analytical Models	Max. weighing capacity	Readability	Linearity	Weighing pan dimensions
APX-60	60 g	0.1 mg	±0.2 mg	Ø 3" (76 mm)
APX-100	100 g	0.1 mg	±0.2 mg	Ø 3" (76 mm)
APX-200	200 g	0.1 mg	±0.2 mg	Ø 3" (76 mm)
Toploading Precision Models				
APX-153	150 g	0.001 g	±0.002 g	Ø 4.5" (114 mm)
APX-203	200 g	0.001 g	±0.002 g	Ø 4.5" (114 mm)
APX-402	400 g	0.01 g	±0.02 g	Ø 4.5" (114 mm)
APX-602	600 g	0.01 g	±0.02 g	Ø 4.5" (114 mm)
APX-1502	1500 g	0.01 g	±0.02 g	Ø 4.5" (114 mm)
APX-2001	2000 g	0.1 g	±0.2 g	5" x 7" (127x178 mm)
APX-4001	4000 g	0.1 g	±0.2 g	5" x 7" (127x178 mm)
APX-6001	6000 g	0.1 g	±0.2 g	5" x 7" (127x178 mm)

General Specifications:

Dimensions (WxDxH), analytical models: 12.7" x 8.8" x 12.25" (322 x 224 x 311 mm)

Dimensions (WxDxH), toploading

precision models:

12.7" x 8.8" x 3.0" (322 x 224 x 76 mm)

Net weight, analytical models:

9 lbs. (4.1 kg) (typical)

Net weight, toploading precision models:

4 lbs. (1.8 kg) (typical)

Power requirements:

15 VDC@100 mA with AC adapter,
center pin (-).

Operating temperature:

15° to 40°C (59° to 104°F)

Storage/transport temperature:

-10° to +40°C (14° to 104°F)

Humidity:

80% for temperatures to 31°C, decreasing
linearly to 50% relative humidity at 40°C

Altitude:

3000 m

Main supply voltage fluctuations must not exceed ±10% of nominal supply voltage. This equipment is suitable for continuous operation with an AC adapter.

Caution:

Use only the AC adapter that is supplied with the balance. Contact the Denver Instrument Company if you need a replacement.

Intended Use

APEX Series balances from Denver Instrument offer precision weighing capacities from 60 to 6000 grams. These balances meet the highest requirements for accuracy and reliability of weighing results with versatile features such as:

- Efficient filtering-out of unfavorable ambient conditions (e.g., vibration or drafts)
- Stable, reproducible results
- Ultra-fast response times
- Easy operation
- Choice of weight units
- Rugged, durable construction suitable for use in training, laboratory and general industrial applications
- Toggling between four weight units, including one user-selected unit, or piece count
- Counting
- Weighing in percent

Warning and Safety Information

Read this operating manual carefully before operating your precision balance. Connect only Denver Instrument accessories and options, as these are designed for optimal performance with your balance.

Warning:

Make sure that the voltage rating printed on the AC adapter is identical to your local line voltage.

Do not use this balance in a hazardous area/location.

The only way to turn power off completely is by disconnecting the AC adapter from the balance.

Protect the AC adapter from contact with liquids.

This product is intended for indoor use.

Warning:

This unit contains no user-serviceable parts. Do not open the balance housing, as this voids the manufacturer's warranty.

Getting Started

Thank you for choosing one of our precision instruments. Your balance is designed and manufactured to the most rigorous standards in order to give you years of dependable service. Please check the contents of the shipping carton as soon as you unpack your balance.

The equipment supplied includes the following components:

- Stainless-steel floor pan (analytical models only)
- Draft shield disk (analytical models only)
- Pan support
- Weighing pan
- AC adapter
- Operating instructions

Next, follow the instructions for installing your balance. To get the most out of its many features, read the operating instructions carefully. They contain step-by-step procedures, examples, and other important information.

Warranty

In the unlikely event that your balance should require servicing within the 3-year warranty period, please contact the office listed below for your country, and state the model number, serial number and the problem:

Federal Republic of Germany:
DENVER INSTRUMENT GmbH
Robert-Bosch-Breite 10
37079 Goettingen
Germany

Telephone: +49 (0) 551/209 773-0

Telefax: +49 (0) 551/209 773-9

Internet:

<http://www.DenverInstrument.com>

E-mail: info@DenverInstrument.de

Unpacking the Balance

After unpacking the balance, check it immediately for any visible damage as a result of rough handling during shipment.

If you see any sign of damage, proceed as directed in the chapter entitled "Care and Maintenance," under "Safety Inspection."

Installation Instructions

Your balance is designed to provide reliable weighing results under normal ambient conditions. When choosing a location to set up your new balance, observe the following so that you will be able to work with added speed and accuracy:

- Set up the balance on an even, stable surface
- Avoid exposing the balance to extreme heat radiation from heaters or direct sunlight
- Do not use the balance at temperatures exceeding the allowable operating temperature range (15°C/59°F to 40°C/104°F)
- Protect the balance from direct exposure to drafts
- Protect the balance from aggressive chemical vapors
- Avoid strong magnetic fields generated by other devices
- Avoid exposing the balance to extreme vibrations during weighing

Preparing the Weighing Chamber

To avoid damage to your precision balance during shipping, the weighing chamber components are packed separately. Toploading models come with a pan support and a weighing pan; analytical models have a separate floor pan and draft shield in addition to the pan support and weighing pan. First install the stainless-steel floor pan, followed by the draft shield (analytical models). Center the pan support over the pan adapter and then position the weighing pan on the pan support.

Connecting the Balance to AC Power



After conditioning the balance, simply plug the AC adapter into the rear of the balance and connect the adapter to an appropriate power outlet. "◇" is shown in the upper left hand corner of the display, indicating that the balance is connected to power and currently in standby mode. Press the I/O key to initialize balance.


Leveling the Balance

Level the balance when first installing and when moving to a new location. Adjust leveling feet until the bubble is centered in the level indicator. Then make sure all feet are touching the countertop. The number of leveling feet varies per model: analytic and toploading models with round pans have 2 feet, toploading models with square pans have 4 feet. In most cases, this requires several adjustment attempts.

Warmup Time


To deliver exact results, the balance must warm up for at least 30 minutes after initial connection to AC power or after a relatively long power outage. Only after this time will the balance have reached the required operating temperature.


Installing the Parts on Balances with a Round Glass Draft Shield

- Attach the contact clip to the top part of the housing.

- Place the draft shield disk on the balance and turn until it audibly clicks into place.
- Install the following components on the balance in the order given:
 - Pan support
 - Weighing pan
 - Glass cylinder
 - Draft shield cover

Operating the Balance

Taring the Balance

1. Press I/O to switch the instrument on. (Allow the balance to warm up to room temperature, a minimum of 30 minutes).
2. A brief self-test is performed, followed by " - - " on the display, indicating that the balance is being zeroed.
3. Place a container on the weighing pan.
The balance registers the weight of the container.
4. Press the Zero key.
5. The display shows " - - " to indicate that the balance is being zeroed.
6. The balance readout shows 0.0000 grams (or other selected weight unit) after successful taring.
7. When the readout is stable, the balance displays the unit of measurement using the following abbreviations:

g	grams	oz	ounces	ozt	Troy ounces
GN	grains	ct	carats	tl	taels
dwt	pennyweights	lb	pounds	pcs	pieces

(in Counting mode)

Toggling between Weight Units

1. You can toggle the weight readout during weighing to show the current weight in grams, pennyweights, ounces, pieces (counting mode) or a user-selected weight unit by pressing the arrow (scroll) key.

Selecting a Weight Unit

1. Press and hold the Print/Menu key for three seconds. The display shows UNIT.
2. Press the arrow (scroll) key repeatedly to scroll through the options (listed above); the display shows cArAt, Lb, tr OZ, tAEL and grAln in turn.
3. When the desired weight unit is displayed, press the ZERO key to return to the weighing mode.

Calibration/Adjustment

“Calibration” technically means to determine the difference between the weighing instrument readout and the actual weight on the balance to determine the accuracy. Adjustment means to bring a weighing instrument to the level of accuracy required for its use.

Always calibrate your balance after it has been disconnected from power, leveled, or moved to a new location. Calibration can be performed only when there is no weight on the balance, the balance is tared and the internal signal is stable. If any of these conditions is not met, an error message is displayed on the screen. The weight required for calibration or adjustment is prompted on the display. (See “Accessories” for calibration weights.)

Calibrating the balance:

1. Press I/O to switch the instrument on. (Allow the balance to warm up to room temperature, approximately 30 minutes).
2. A brief self-test is performed, followed by “- -” indicating that the balance is being zeroed.
3. Remove all items from the pan and press the Zero key to tare balance.
4. Press and hold the Print/Menu key for three seconds or until “Unit” is displayed.
5. Press the Print/Menu key again; “CAL” is displayed.
6. Press the arrow key to activate the calibration mode.
7. The required calibration weight is displayed (e.g., +200.0000) in grams.
8. Place the prompted calibration weight on the pan.
9. When the readout has stabilized, a beep-tone is emitted, the “- -” sign flashes on the display, and the readout returns to the display of the calibration weight.
10. Remove the calibration weight; you can now begin weighing.

Counting

With the Counting program you can determine a number of parts, each having approximately the same weight. The weight of 10, 20, 50 or 100 pieces is stored as the reference sample quantity first, and the next load on the balance is then calculated as a number of pieces having the same weight, based on this reference value.

Available reference sample quantities:

- 10 pieces
- 20 pieces
- 50 pieces
- 100 pieces

Factory default:

Reference sample quantity: 10 pieces

Activating the counting application and setting the reference sample quantity:

1. Press I/O to switch the instrument on.
2. A brief self-test is performed, followed by “- -” indicating that the balance is being zeroed.
3. Press and hold the Print/Menu key for three seconds or until “Unit” is displayed.
4. Press the Print/Menu key twice; “Count” is displayed.
5. Press the arrow key to activate the counting mode.
6. Select the reference sample quantity by pressing the arrow key repeatedly until the desired quantity is displayed. The display shows: 10, 20, 50, 100, and P100.0 in turn.
7. Press the Zero key to confirm.
8. Remove all items from the pan and press the Zero key to tare the balance.
9. Place the reference quantity (e.g., 10 pieces) on the balance and wait for the readout to stabilize (the weight unit is displayed).

10. Press the arrow key. The balance stores the current value as the reference weight and the display shows the current piece count.
11. Add uncounted parts. A total number of pieces is displayed.
12. Zero the balance to change the piece count to zero.
13. Press the arrow key (repeatedly, if necessary) to display the total weight of pieces in grams, penny-weights, ounces, or in the user-selected unit.

The reference weight is reset when you toggle from total number of parts counted to total weight.

Changing the reference weight:

1. Press the arrow key repeatedly until the weight of the pieces is displayed.
2. Remove all items from the pan and add a number of pieces corresponding to the reference sample quantity.
3. Press the arrow key repeatedly until “- -” is displayed, indicating that the balance is storing the weight as the new reference weight. The total number of pieces is now displayed.
4. Add uncounted parts. The total number of pieces is displayed.

When scrolling through weight units, the counting mode is displayed, but the reference sample quantity is reset only if there is weight on the pan.

Weighing in Percent

With the Weighing in Percent program, you can obtain weight readouts in a percentage proportionate to a reference weight. A weight on the balance is stored as the reference weight (100%) first, and the next load on the balance is then measured as a percentage of this weight.

Setting the reference weight:

1. Press I/O to switch the instrument on.
2. A brief self-test is performed, followed by "- -" indicating that the balance is being zeroed.
3. Press and hold the Print/Menu key for three seconds or until "Unit" is displayed.
4. Press the Print/Menu key twice; "Count" is displayed.
5. Press the arrow key to activate the counting mode.
6. Press the arrow key repeatedly until P100.0 is displayed. The display shows: 10, 20, 50, 100, and P100.0 in turn.
7. Press the Zero key to confirm.
8. Remove all items from the pan and press the Zero key to tare the balance.
9. Place the reference weight on the balance and wait for the readout to stabilize (indicated by display of the weight unit).
10. Press the arrow key repeatedly until the display shows 100.00%.
11. Add or remove parts; a percentage value is displayed.
12. Press the arrow key (repeatedly, if necessary) to display the total weight in grams, pennyweights, ounces, or in the user-selected unit.

The reference weight is reset when you toggle from percentage ("%") to total weight.

Changing the reference weight:

1. Press the arrow key repeatedly until the weight of the pieces is displayed.
2. Remove all items from the pan and add the new total ("100%") weight.
3. Press the arrow key repeatedly until "- -" is displayed, indicating that the balance is storing the weight as the new reference weight. 100% is now displayed.
4. Add or remove parts; a percentage value is displayed.

When scrolling through units, the percent mode is displayed and reset only when there is weight on the pan.

Configuring the Balance

Several advanced balance operations can be changed through menu codes. Select the desired codes from the list below. Then follow the steps to change the menu codes.

Special key functions for menu code operation:

Zero- To enter/leave menu code settings; increases a number by one with each press

Print- Moves to the next of three numbers of a code

To Change Menu Codes:

1. Press power key to turn balance off
2. Press power key to turn the balance back on. While all segments are displayed, press zero key once. After self-test, screen will have 1 on the left side of display
3. Press zero key to increase the number on display to desired value
4. Press print for APX to move to the 2nd number of the code
5. Press zero key to increase the second number on display to desired value
6. Press print for APX to move to the 3rd number of the code (when you move to the third number, the previously set menu code will appear)
7. Press zero key to increase the third number on display to desired value
8. Hold down print for APX until "o" appears on the APX after the set code. This confirms and stores the new code.
9. Press zero for more than 2 seconds to store the new menu code setting.

Menu Codes:

Filter Selection

- 1 1 1 Very stable conditions
- 1 1 2 Stable conditions (default)
- 1 1 3 Unstable conditions
- 1 1 4 Very unstable conditions

Baud rate

- 5 1 1 150 baud
- 5 1 2 300 baud
- 5 1 3 600 baud
- 5 1 4 1200 baud
- 5 1 5 2400 baud
- 5 1 6 4800 baud
- 5 1 7 9600 baud (default)

Parity

- 5 2 1 Mark
- 5 2 2 Space (default)
- 5 2 3 Odd
- 5 2 4 Even

Stop Bits

- 5 3 1 1 stop bit
- 5 3 2 2 stop bits


Handshake Mode


- 5 4 1 Software
- 5 4 2 Hardware, 2 characters after CTS
- 5 4 3 Hardware, 1 character after CTS

Print Configuration

- 6 1 1 Manual with print key without stability
- 6 1 2 Manual with print key after stability (default)
- 6 1 3 Manual with print key at stability
- 6 1 4 Auto print without stability

Device Parameters (Overview)

	Function in Weighing Mode	Function in Menu Mode
	Zeroes the balance	Exits current weighing mode, saves (new) settings
	Pressed briefly: sends printer string to optional serial port. Pressed for 3 seconds: opens the operating menu	Advances to next menu line
	Toggles the display mode through all options in series	Scrolls through items on current menu level

Notes:

1. Texts in quotes above describe the functions activated by each key (e.g., "return to weighing mode")
2. Texts not in quotes are display messages (e.g., count or tael) or describe user actions (e.g., "place weight")
3. P100.0 denotes percentage weighing

Functions of each menu line:

- unit define the weight unit available in addition to the default units (grams, ounces, pennyweights, piece count)
- cal calibrate the weighing instrument to the applied reference weight
- count select the reference sample quantity and activate the counting mode, or select the reference percentage (P100.0) and activate the weighing in percent mode

Troubleshooting

Display:	Cause:	Solution:
No segments appear on the display.	Display switched off by pressing I/O key.	Press the I/O key.
No \diamond displayed	Power cord not connected.	Connect cord.
	No power to outlet or incorrect voltage.	Check power supply. Check power supply and voltage switch.
----	The load exceeds the balance capacity.	Unload the balance or remove obstruction.
	Display capacity is exceeded.	Decrease weight on balance.
E 22	The weighing pan is not positioned correctly.	Position the weighing pan correctly or remove obstruction
	Weight is too light or there is no sample on balance.	Increase reference quantity.
E 01	Display capacity exceeded, value to be output cannot be displayed.	Decrease weight on the balance.
E 02	Calibration condition not met.	Calibrate only when zero is displayed.
	Balance not tared.	Press ZERO key.
	Balance is loaded.	Unload balance.
E 11	Value input is not allowed for second tare memory.	Press ZERO key.
E 30	Interface port for printer output is blocked.	Contact Denver Instrument Company.
Maximum weighing range is less than indicated under "Specifications"	The balance was switched on without the weighing pan in place.	Place the weighing pan on the balance and press the I/O key to switch balance on.
The weight readout is obviously incorrect.	Balance has not been calibrated.	Calibrate the balance.
	Balance was not tared.	Tare balance.

If other errors occur, contact your Denver Instrument Service Representative.

RS-232C Interface

Data Input Format

You can connect a computer or printer to your balance to send commands via the balance interface port to control balance functions and applications.

A control command can have up to 4 characters. Each character must be transmitted according to the interface port settings for data transmission.

Interface port settings:

Baud	Data bits	Parity	Stop bit
9600	7	Space	1

Software handshake: =	None, or Xon/Xoff
Line output interval:	50 ms
Character output interval:	50 ms

Format for control commands:

Format: Esc ! CR LF

Esc: Escape	CR: Carriage Return (optional)
!: Command character	LF: Line Feed (optional)

Command character:	!	Meaning
	K	Very stable conditions
	L	Stable conditions
	M	Unstable conditions
	N	Very unstable conditions
	O	Disable keypad
	P	Print display
	R	Enable keypad
	S	Reset balance
	T	Tare and zero (combined)
	U	Tare (tare only)
	V	Zero
	W	External calibration/adjustment

Hardware Handshake

With a 4-conductor interface, 1 or 2 characters can be transmitted after CTS.

These connections must be made when the balance is connected through the RS-232C port.

Pin Assignment Chart for Data Interface

APEX Series		Standard RS-232
Balance 9-pin port		9-pin connector
RxD	2	3 TxD
TxD	3	2 RxD
DTR	4	4 DTR
Signal ground	5	5 Signal ground
CTS	8	8 CTS

APEX Series		Standard RS-232
Balance 9-pin port		25-pin connector
RxD	2	2 TxD
TxD	3	3 RxD
DTR	4	20 DTR
Signal ground	5	7 Signal ground
CTS	8	5 CTS

Care and Maintenance

Service

Regular service by a Denver Instrument technician will extend the life of your balance and ensure its continued weighing accuracy.

Repairs

Repair work must be performed by qualified factory-trained personnel only.

Note:

This unit contains no user serviceable parts.

All replacement parts should be obtained from the manufacturer. Please refer to the inside front cover of this manual for the phone number of your sales and service representative.

Warning:

Never lift the balance by the draft shield as this may cause damage to internal mechanisms.

Always lift and transport the balance by its base.

Cleaning

Caution: Disconnect the balance AC adapter from power source prior to cleaning. Make sure that no liquids enter the balance housing.

Do not use aggressive or abrasive cleaning agents, such as cleansers.

A mild detergent is recommended.

Disassemble the weighing chamber (see page 9) and clean the floor pan, draft shield, pan support and weighing pan separately, then reassemble.

Use a piece of cloth to clean the balance. After cleaning, wipe the balance down with a dry, soft cloth. Recalibration of the balance is recommended after cleaning.

Safety Inspection

Safe operation of the balance with the AC adapter is no longer ensured if there is visible damage to the AC adapter or cord, the AC adapter no longer functions properly or the AC adapter has been stored for a relatively long period under unfavorable conditions.

Accessories

Draft shield disk

Dot matrix printer

9-pin cable

25-pin cable

Calibration weight – 50 g

Calibration weight – 50 g

Calibration weight – 100 g

Calibration weight – 200 g

Calibration weight – 500 g

Calibration weight – 1000 g

Calibration weight – 2000 g

902191.1

Information available on request

Information available on request

Information available on request

Information available on request

Information available on request

Information available on request

Information available on request

Information available on request

Information available on request

Information available on request

**Declaration of
Conformity**

**to Council
Directives
89/336/EEC
and 73/23/EEC**

**DENVER
INSTRUMENT**

**The electronic precision weighing
instrument of the series**

APX-....

meets the requirements of the test standards listed below, in conjunction with the associated power supplies, auxiliary peripheral devices and installation equipment listed in Annex A2 (see Annex A1 for a technical description and variants).

1. Electromagnetic Compatibility

1.1 Source for 89/336/EEC: EC Official Journal, No. 2001/C105/03

EN 61326-1 Electrical equipment for measurement, control and laboratory use
EMC requirements
Part 1: General requirements

Limitation of emissions: Residential areas, Class B
Defined immunity to interference: For use in controlled EM environments only

2. Safety of Electrical Equipment

2.1 Source for 73/23/EEC: EC Official Journal, No. 2001/C106/03

EN 61010 Safety requirements for electrical equipment for measurement, control and laboratory use
Part 1: General requirements

EN 60950 Safety of information technology equipment

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