



IC System Operational and Performance Qualification User's Guide

© 2000 Dionex Corporation

Document No. 031726 - 02
November 2000

© 2000, Dionex Corporation
All rights reserved worldwide.
Printed in the United States of America

This publication is protected by federal copyright law. No part of this publication may be copied or distributed, transmitted, transcribed, stored in a retrieval system, or transmitted into any human or computer language, in any form or by any means, electronic, mechanical, magnetic, manual, or otherwise, or disclosed to third parties without the express written permission of Dionex Corporation, 1228 Titan Way, Sunnyvale, California 94088-3606 USA.

DISCLAIMER OF WARRANTY AND LIMITED WARRANTY

THIS PUBLICATION IS PROVIDED "AS IS" WITHOUT WARRANTY OF ANY KIND. DIONEX CORPORATION DOES NOT WARRANT, GUARANTEE, OR MAKE ANY EXPRESS OR IMPLIED REPRESENTATIONS REGARDING THE USE, OR THE RESULTS OF THE USE, OF THIS PUBLICATION IN TERMS OF CORRECTNESS, ACCURACY, RELIABILITY, CURRENCY, OR OTHERWISE. FURTHER, DIONEX CORPORATION RESERVES THE RIGHT TO REVISE THIS PUBLICATION AND TO MAKE CHANGES FROM TIME TO TIME IN THE CONTENT HEREIN WITHOUT OBLIGATION OF DIONEX CORPORATION TO NOTIFY ANY PERSON OR ORGANIZATION OF SUCH REVISION OR CHANGES.

PRINTING HISTORY

Revision 01, August 2000
Revision 02, November 2000

Table of Contents

1	Introduction	4
1.1	Operational Qualification (OQ)	5
1.2	Performance Qualification (PQ)	5
1.3	OQ/PQ Specifications	6
2	Qualification Preparation	7
2.1	Materials in Kit	7
2.2	Materials Needed	8
2.3	Preparation of Standards	8
2.3.1	Caffeine Standard Preparation	8
2.3.2	Nitrate Standard Preparation	9
2.4	Preparation of Eluents	9
2.4.1	Preparation of Eluents for use with Conductivity Detection	9
2.4.2	Preparation of Eluents for use with Absorbance Detection	10
2.5	IC System Preparation	10
2.6	OQ/PQ Sequence Installation	13
2.6.1	Installing Sequences from the PeakNet 6.20 CD	13
2.6.2	Installing Sequences from the OQ/PQ CD (P/N 057601)	13
2.7	OQ/PQ Report Set Up	17
3	Performing the Qualification	18
3.1	When to Perform Operational Qualification	18
3.2	When to Perform Performance Qualification	18
3.3	OQ/PQ Test Order	18
3.4	What to Do if a Test Fails	18
3.5	Starting the Sequences	18
3.5.1	Wavelength Accuracy	19
3.5.2	Performing the Main Qualification	21
3.5.3	Gradient Accuracy	22
3.5.4	DC Amperometry Noise and Drift	23
3.5.5	Integrated Amperometry Noise and Drift	24
3.5.6	EG40 Current Qualification	24

3.6	OQ/PQ Completion	26
3.6.1	Customer Review	26
3.6.2	Returning the System to the Original Configuration	26
4	Explanation of Test Procedures	27
4.1	Wavelength Accuracy of Absorbance Detectors	27
4.1.1	Theory	27
4.1.2	Performing and Evaluating the Test	27
4.2	Baseline Noise and Drift	27
4.2.1	Theory	27
4.2.2	Performing the Test and Evaluating the Results	27
4.3	Injector Precision	28
4.3.1	Theory	28
4.3.2	Performing and Evaluating the Test	28
4.4	Detector Linearity	28
4.4.1	Theory	28
4.4.2	Performing the Test and Evaluating the Results	28
4.5	Injection Volume Linearity	29
4.5.1	Theory	29
4.5.2	Performing the Test and Evaluating the Results	29
4.6	Sample Carry-over	29
4.6.1	Theory	29
4.6.2	Performing the Test and Evaluating the Results	29
4.7	Sample Prep Precision	30
4.7.1	Theory	30
4.7.2	Performing the Test and Evaluating the Results	30
4.8	Gradient Accuracy	30
4.8.1	Theory	30
4.8.2	Performing the Test and Evaluating the Results	30

4.9	DC Amperometry Noise and Drift	31
4.9.1	Theory	31
4.9.2	Performing the Test and Evaluating the Results	31
4.10	Integrated Amperometry Noise and Drift	31
4.10.1	Theory	31
4.10.2	Performing the Test and Evaluating the Results	31
4.11	EG40 Current Qualification	31
4.11.1	Theory	31
4.11.2	Performing the Test and Evaluating the Results	32
5	Troubleshooting	33
5.1	General Hints	33
5.2	Failure of Individual Tests	34
5.2.1	Absorbance Detector	34
5.2.2	Conductivity Detector Tests	35
5.2.3	Autosampler	36
5.2.4	Pump	37
5.2.5	DC and Integrated Amperometry Detector Tests	37
5.2.6	EG40	38
6	Modifying an EG40 Validation Test Cell	39
7	Example Reports	40

1 Introduction

Validation is becoming increasingly important to analytical laboratories. Documented evidence must be provided to demonstrate the integrity of data collected and validate the results obtained on laboratory instrumentation.

The following institutes set and control industry instrumentation standards:

NIST. The National Institute of Standards and Technology is an administrative body of the US Department of Commerce. Together with the industry, the NIST develops new technologies and sets technological standards countrywide.

ASTM. The American Society for Testing and Materials sets technical standards and rules for the industry worldwide.

These standards are written in broad terms to make them as widely applicable as possible. All stipulate general requirements specifying instruments must be fit for purpose, technically checked in defined periods, calibrated and make sure that test results are documented according to international quality standards. The OQ/PQ procedure used to qualify Dionex Ion Chromatography Systems meets the requirements established by NIST and ASTM.

This QQ/PQ procedure provides qualification testing for **DX-500, DX-600, DX-320 and BioLC** systems.

OQ/PQ results can be influenced by unusual laboratory conditions. For example, direct sunlight will cause thermal irregularities that will affect baseline drift on conductivity detectors. In these cases it may happen that the manufacturer's specifications are not met, and the specifications may have to be adjusted accordingly.

All OQ/PQ activities must be formally recorded, all instruments must be clearly identified, and test reports and supporting documentation must include the following information:

- Module name, type and serial number
- Software and firmware version
- Test data
- Service interval
- Signature of the Dionex Service Engineer/Qualification Executor who performed the OQ/PQ
- Signature of the customer who reviewed and accepted the OQ/PQ results

1.1 Operational Qualification (OQ)

The purpose of OQ is to test a new instrument according to the manufacturer's specifications within a specific environment. Therefore, the manufacturer must define the test conditions and the test specifications, both which must be observed.

OQ should be performed after initial system installation or when an additional new module is added to an existing system. It is necessary that OQ be done according to the instructions outlined in this manual and all tests should be performed by a qualified Dionex Service Engineer.

1.2 Performance Qualification (PQ)

The purpose of PQ is to determine if a system that has been subjected to daily routine operation functions within the Dionex recommended limits of acceptable system performance.

Since PQ is performed on a system that is no longer considered new, the specifications for PQ are less restrictive than those for OQ, even though the same test procedures are used for both. PQ results must be either within the limits recommended by Dionex or be adapted to the limits required by the customer. If customer-specified limits are used, those limits must not be tighter than the Dionex specification.

PQ should be performed at regular intervals after the initial installation and OQ. Dionex recommends performing PQ every six months. It is necessary that PQ be done according to the instructions outlined in this manual and all tests should be performed by a qualified Dionex Service Engineer.

1.3 OQ/PQ Specifications

Instrument	Test	Limits	
		OQ	PQ
AD25	Baseline Noise	30 μ AU	40 μ AU
	Drift	200 μ AU/hr	200 μ AU/hr
	Wavelength Accuracy	+/- 2.00 nm	+/- 2.00 nm
	Linearity	$R^2 \geq 0.999$ RSD $\leq 5\%$	$R^2 \geq 0.999$ RSD $\leq 5\%$
AD20	Baseline Noise	40 μ AU	50 μ AU
	Drift	200 μ AU/hr	200 μ AU/hr
	Wavelength Accuracy	+/- 3.00 nm	+/- 3.00 nm
	Linearity	$R^2 \geq 0.999$ RSD $\leq 5.0\%$	$R^2 \geq 0.999$ RSD $\leq 5.0\%$
Conductivity CD20/ED40, CD25/ED50, IC20/IC25	Baseline Noise	1.0 nS	2 nS
	Drift	10 nS/hr	20 nS/hr
	Linearity	$R^2 \geq 0.999$ RSD $\leq 5.0\%$	$R^2 \geq 0.999$ RSD $\leq 5.0\%$
Integrated Amperometry	Baseline Noise	3 pC	4 pC
	Drift	3 pC/hr	5 pC/hr
DC Amperometry	Baseline Noise	3 pA	4 pA
	Drift	4 pA/hr	8 pA/hr
PDA -100	Baseline Noise	0.04 mAU	0.10 mAU
	Drift	0.5 mAU/h	1.0 mAU/h
	Wavelength Accuracy	+/- 2.00 nm	+/- 2.00 nm
	Linearity	$R^2 \geq 0.999$ RSD $\leq 5.0\%$	$R^2 \geq 0.999$ RSD $\leq 5.0\%$
AS50	Injector Precision	RSD $\leq 0.8\%$	RSD $\leq 1.0\%$
	Injection Volume Linearity	$R^2 \geq 0.999$ RSD $\leq 1.0\%$	$R^2 \geq 0.999$ RSD $\leq 1.0\%$
	Carry Over	< 0.1 %	< 0.1 %
	Sample Prep Precision	RSD $\leq 1.5\%$	RSD $\leq 2.0\%$
GP40/GP50 Standard Bore, GS50	Gradient Accuracy	$\leq 1.5\%$	$\leq 2.0\%$
	Step Noise (Ripple)	$\leq 0.2\%$	$\leq 0.2\%$
GP40/GP50 Microbore	Gradient Accuracy	$\leq 2.0\%$	$\leq 2.0\%$
	Step Noise (Ripple)	$\leq 0.5\%$	$\leq 0.5\%$

Note: The instrument's specifications can only be met if the conditions described in the accompanying manual are strictly observed.

2 Qualification Preparation

2.1 Materials in Kit

Before starting the qualification procedure confirm that all the materials needed are present. If any items from these tables are not present, do not continue with the validation process until these items are available.

Part Description	Part Number	Quantity
OQ/PQ User's Manual Version 1.0	031726	1
Tubing, 0.010" (0.25 mm) ID black PEEK	042690	24", 81 cm
Tubing, 0.003" (0.075 mm) ID yellow PEEK	049715	50", 127 cm
25 µL Sample Loop	042857	1
100 µL Sample Loop	042951	1
Fittings	043275	6
Ferrules	043276	6
Validation stickers	050837	6
Certificate of Functional Performance	031423	6
1000 ppm Caffeine Solution	056822	1
1000 ppm Nitrate Solution	056497	1
Backpressure Coil, 0.005" (0.125 mm) ID PEEK x 512" (13 m)	2251.5999	1
ED/CD Validation Test Cell*	049928	1
EG40 Validation Test Cell*	054970	1
AD20 Holmium Oxide Test Cell*	046172	1
OQ/PQ Validation Software for CM/PN 6	057601	1
Test Binder	035863	1

* Included in the PeakNet 6 IC OQ/PQ Kit with Test Cells (Dionex P/N 057599) only. If using the PeakNet 6 IC OQ/PQ Basic/Refill Kit (Dionex P/N 057608), the Dionex Service Engineer must provide test cells.

2.2 Materials Needed

The following items will be needed:

Description	Quantity
ASTM Type 1 (or better) deionized water (18 M Ω , filtered)	2 L
Reagent Grade or better Potassium Chloride for OQ/PQ tests using a Conductivity Detector	1 g
HPLC grade Acetone for OQ/PQ Tests using an Absorbance or Photodiode Array Detector	3 mL
Voltmeter with milliamp (mA) setting. Must be capable of displaying 0.00 to 220 mA.	1
PeakNet [™] Version 6.01 or higher	1

2.3 Preparation of Standards

Two standards are provided in the OQ/PQ kit. Use the caffeine standard to qualify systems with an AD20, AD25 or PDA-100 detector. Use the nitrate standard to qualify systems with an ED40/ED50, CD20/CD25 or IC20/IC25 Conductivity detector. Only one set of standards needs to be prepared depending on the type of detector used in the qualification.

2.3.1 Caffeine Standard Preparation

Prepare caffeine standards using the Dionex 1000 ppm caffeine standard (P/N 056822). Pipette the amount of 1000 ppm caffeine standard (listed in the table below) into five individual 100 mL volumetric flasks. Fill each flask to the mark with ASTM Type 1 (or better) deionized water (18 M Ω , filtered).

Concentration (ppm)	Amount of 1000 ppm Caffeine Standard	Final Volume (mL)
15	1.5 mL	100
40	4.0 mL	100
60	6.0 mL	100
80	8.0 mL	100
100	10.0 mL	100

2.3.2 Nitrate Standard Preparation

Prepare nitrate standards using the Dionex 1000 ppm nitrate standard (P/N 056497). Pipette the amount of 1000 ppm nitrate standard listed in the table below into five individual 100 mL volumetric flasks. Fill each flask to the mark with ASTM Type 1 (or better) deionized water (18 MΩ, filtered).

Concentration (ppm)	Amount of 1000 ppm Caffeine Standard	Final Volume (mL)
5	0.5 mL	100
10	1.0 mL	100
25	2.5 mL	100
50	5.0 mL	100
100	10.0 mL	100

2.4 Preparation of Eluents

All eluents must be prepared using ASTM Type 1, 18 MΩ or better, filtered, deionized water.

2.4.1 Preparation of Eluents for use with Conductivity Detection

These eluents are used for qualification with an ED40/ED50, CD20/CD25 or IC20/25 Conductivity Detector.

Prepare the following eluents and connect them to lines A, B, C and D of the pump. All eluents must be degassed. Eluent bottles must be pressurized when using a GP40 or IP20 pump. Eluent lines A and C can be inserted into one eluent bottle containing 1 L of water. Eluent lines B and D can be inserted into one bottle containing 1 L of 1 mM KCl.

Eluent	Description
A	Water
B	1.0 mM KCl
C	Water
D	1.0 mM KCl

Note: Eluents B, C and D are needed for the Gradient Pump Accuracy Test only and should be prepared only if the pump configured has gradient capabilities. Do not prepare these eluents if the pump configured is an IP20/IP25, IC20/IC25, or IS25 Isocratic pump.

1 mM KCl Preparation:

Prepare one liter of 10 mM KCl by weighing 0.746 grams of KCl. Place the KCl in a 1L volumetric flask. Fill flask to the mark with ASTM Type 1 (or better) deionized water (18 MΩ, filtered) water.

Pipette 100 mLs of the 10 mM KCl solution into another 1L volumetric flask, filling flask to the mark with ASTM Type 1 (or better) deionized water (18 MΩ, filtered) water.

2.4.2 Preparation of Eluents for use with Absorbance Detection

These eluents are for used for qualification with an AD25 or AD20 Absorbance Detector or the PDA-100 Photodiode Array Detector.

Prepare the following eluents and connect them to lines A, B, C and D of the pump. All eluents must be degassed. Eluent bottles must be pressurized when using a GP40 or IP20 pump. Eluent lines A and C can be inserted into one eluent bottle containing 1 L of water. Eluent lines B and D can be inserted into one bottle containing 1 L of 0.3% Acetone in Water.

Eluent	Description
A	Water
B	0.3% Acetone in Water
C	Water
D	0.3% Acetone in Water

Note: Eluents B, C and D are needed for the Gradient Pump Accuracy Test only and should be prepared only if the pump configured has gradient capabilities. Do not prepare these eluents if the pump configured is an IP20/IP25, IC20/IC25, or IS25 Isocratic pump.

0.3% Acetone in Water Preparation:

Pipette 3 mLs of HPLC grade Acetone into a 1L volumetric flask. Fill flask to the mark with ASTM Type 1 (or better) deionized water (18 M Ω , filtered) water.

2.5 IC System Preparation

Note: The AD25 and PDA-100 Detector UV and Visible lamps must be turned on at least two hours prior to performing the qualification. AD20 Detector UV and visible lamps must be turned on for at least 24 hours before performing the qualification.

1. Open Server Configuration. Double click on the detector/s to be qualified to view the properties page for that detector. Select the **Signals** page and confirm that the signal factor is 1.000. If it is not, change the factor to 1.000 for the qualification. Exit the properties page and save the configuration.
2. Connect Eluents A, B, C and D to the pump and prime each eluent line individually for approximately 3 minutes each, longer if the vacuum degas option is installed. Eluent bottles must be pressurized when using GP40 and IP20 pumps.
3. Disconnect the column and suppressor from the system and thoroughly rinse all fluid connections with water for approximately 15 minutes.
4. Leaving the gradient mixer in-line, connect the pump line directly to port 2 on the injection valve. If an EG40 is in-line with the IC system, bypass all plumbing to the EG40 and connect the pump directly to the injection valve.

5. Connect one end of the Backpressure Coil (P/N 2251.5999) to port 3 on the injection valve and the other end directly to the inlet of the conductivity or absorbance cell, bypassing all columns and suppressors.
6. Thoroughly rinse the autosampler fluidics with water. This can be accomplished on the AS50 by flushing the sample syringe two or three times and on the AS40 by running two or three rinse vials.
7. For systems using an AS50 Autosampler, install the 100 μ L sample loop on the injection valve. Then, from the front panel of the AS50, go to the **Main Menu**, select **Module Set-Up** and then **Plumbing Configuration**. Change the sample loop size to 100 μ L. For systems using an AS40 Autosampler or manual injector, install the 25 μ L sample loop on the injection valve.
8. Select 100% of Eluent A (water) and a flow rate of 1.00 mL min. Turn the pump flow on and let the system equilibrate for a minimum of 15 minutes.
9. Check that there are no fluctuations in pressure when changing from 'Load' to 'Inject' and vice versa. Failure for the pressure to remain stable when the injection valve actuates could indicate a leak, blockage, or contamination in the system, which must be eliminated before starting the qualification.
10. If using an AS40 autosampler, connect the 'LOAD IN' relay on the rear panel of the AS40 to RLY 1 OUT on the pump.
11. Transfer the standards prepared in Section 2.3 into autosampler vials and place them in the following position/order according to the tables below. To qualify an IC system by manual injection, use the AS40 table for the injection order.

Absorbance or Photodiode Array Detector with AS50	Conductivity Detector with AS50	Vial Tray Position
15 ppm Caffeine Std	5 ppm Nitrate Std	1
40 ppm Caffeine Std	10 ppm Nitrate Std	2
60 ppm Caffeine Std	25 ppm Nitrate Std	3
80 ppm Caffeine Std	50 ppm Nitrate Std	4
100 ppm Caffeine Std	100 ppm Nitrate Std	5
1000 ppm Caffeine Std	1000 ppm Nitrate Std	6
Water	Water	7
Empty Vial*	Empty Vial*	8
Empty Vial*	Empty Vial*	9
Empty Vial*	Empty Vial*	10

* Needed only if qualifying Sample Prep Option. AS50 10 mL vials are required for Sample Prep Qualification.

Absorbance or Photodiode Array Detector with AS40	Conductivity Detector with AS40	Sample Order
40 ppm Caffeine Std	50 ppm Nitrate Std	Vial 1
40 ppm Caffeine Std	50 ppm Nitrate Std	Vial 2
40 ppm Caffeine Std	50 ppm Nitrate Std	Vial 3
40 ppm Caffeine Std	50 ppm Nitrate Std	Vial 4
40 ppm Caffeine Std	50 ppm Nitrate Std	Vial 5
40 ppm Caffeine Std	50 ppm Nitrate Std	Vial 6
15 ppm Caffeine Std	5 ppm Nitrate Std	Vial 7
40 ppm Caffeine Std	10 ppm Nitrate Std	Vial 8
60 ppm Caffeine Std	25 ppm Nitrate Std	Vial 9
80 ppm Caffeine Std	50 ppm Nitrate Std	Vial 10
100 ppm Caffeine Std	100 ppm Nitrate Std	Vial 11
1000 ppm Caffeine Std	1000 ppm Nitrate Std	Vial 12
Water	Water	Vial 13

2.6 OQ/PQ Sequence Installation

2.6.1 Installing Sequences from the PeakNet 6.20 CD

1. Create a directory in the local datasource called **Validation**.
2. Insert the PeakNet 6.20 CD into the computer's CD drive
3. From the PeakNet Browser, go to **File** and then to **Mount Datasource**, and select **Browse**.
4. From the PeakNet 6.20 CD, select the **IC_OQ_PQ_R02** directory. Select **cm_local.mdb**. The IC_OQ_PQ_R02 datasource will now be available from the PeakNet Browser. Open the datasource, and then open either the IC_OQ or IC_PQ directory, depending on the qualification that will be performed.
5. Proceed to Step 8 in Section 2.6.2.

2.6.2 Installing Sequences from the OQ/PQ CD (P/N 057601)

The OQ/PQ CDROM (Dionex P/N 057601) is included in the PeakNet 6 IC OQ/PQ Kit with Test Cells (Dionex P/N 057599) and the PeakNet IC OQ/PQ Basic/Refill Kit (Dionex P/N 057608).

1. Create a directory in the local datasource called **Validation**.
2. Insert the OQ/PQ CD into the computer's CD drive.
3. Open Windows Explorer and confirm that the PeakNet cmb file names, file size and modified date/time on the CD matches those in the table below:

File Name	Size	Modified Date	Modified Time
IC_OQ_R02.cmb	1369 KB	11/6/00	4:29 PM
IC_PQ_R02.cmb	1369 KB	11/6/00	4:31 PM

4. From the PeakNet Browser, highlight the local datasource (i.e.Dionex_local). Then select **File, Import/Restore, and Restore...**
5. Open the **IC_OQ.cmb** or **IC_PQ.cmb** file depending on the qualification that will be performed.
6. Select **Other** in the **Destination** box. The local datasource should be displayed. If it is not, use the Browse function to select it. Click on **OK** to install the cmb file into the local datasource.

7. Right click on the newly installed IC_OQ or IC_PQ directory and select **Properties**. Go to the Statistics tab and confirm that the statistics in the **Amount** column match the screen below. Both cmb files will have these statistics:

Access Control		Statistics	
Object of Type	Amount	Disk Space	

Sub-Directories	8		
Sequences	73		
Samples (Total)	605		
Samples (without time stamp)	605		

Raw Data Files			
Chromatograms	0	0 KByte	
3D-Fields	0	0 KByte	
Mass Spectrometry (Spectra)	0	0 KByte	
Audit Trails	0	0 KByte	

Programs (PGM-Files)	276	342 KByte	
Quantification Methods (QNT-Files)	123	309 KByte	
Sample Queries (QRF-Files)	0	0 KByte	
Spectra Libraries (LIB-Files)	0	0 KByte	
Report Definitions (RDF-Files)	2	2439 KByte	
Online Panels (PAN-Files)	0	0 KByte	
Workspaces (WSP-Files)	0	0 KByte	
Signed Results (SOR-Files)	0	0 KByte	
MS Instr. Methods (METH-Files)	0	0 KByte	

History Entries	0		
Detail History Entries	0		

Total occupied disk space:		3090 KByte	

8. Before performing the qualification, sequence templates must be selected and copied from the **IC_OQ** or **IC_PQ** directory into the **Validation** directory. Having the sequences in one location will make them easier to find when starting the batch in Section 3.

9. The OQ/PQ qualification tests are performed in the following order:

- A. Wavelength Accuracy Test for Absorbance Detectors
- B. Detector Noise and Drift Test
- C. Injector Precision Test
- D. Detector Linearity Test
- E. Injector Linearity Test - *AS50 only*
- F. Carry Over Test
- G. Sample Prep Precision Test - *AS50 with Sample Prep Option only*
- H. Gradient Accuracy Test*
- I. DC Amperometry Baseline Noise and Drift*
- J. Integrated Amperometry Baseline Noise and Drift*
- K. EG40 Current Qualification Test*

* Perform this test only if module is configured and/or used in this mode on the timebase to be qualified.

10. Test A is done manually and does not require a sequence template. Sequence templates must be selected for all other tests. Tests B - G are contained in one sequence template. Tests H - K have individual sequence templates.

11. To begin copying sequences, a Detector_Pump sequence template must be selected. The Detector_Pump sequence performs Tests B-G.

Note: Full qualification can only be performed with the **AD20, AD25, PDA-100** or **Conductivity** Detectors. If the system to be qualified only contains an Integrated Amperometry or DC Amperometry detector, this part of the qualification can not be performed. Go on to step 16.

12. Open the **DX320** or **DX500_DX600_BioLC** directory depending on the system to be qualified.

13. Next open the **AS40 and Man Inj** or **AS50** directory, depending on the autosampler or injector configured on the timebase to be qualified. For systems configured with manual or automated injection valves, use the AS40 sequences. For systems configured with an AS50, open the AS50_TC or AS50_CC directory and choose a standard or sample prep (SP) sequence.

Example: If the system to be qualified consists of an AD25, an IP25 and an AS50 with Sample Prep and CC options, you would select the **AD25_CC_SP_Pump** sequence.

Example: If the system to be qualified consists of a CD20, a GP40 and an AS40, you would select the **ED_CD_AS40_Pump** sequence

Example: If the system to be qualified consists of an AD20 and a GP50 with a manual or automated injection valve, you would select the **AD20_AS40_Pump** sequence.

14. Right click on the Detector_Pump sequence template and select **Copy**. Right click on the **Validation** directory and select **Paste**.

15. If a gradient pump is to be qualified, select the appropriate **Grad Accuracy** sequence template and copy it to the **Validation** directory. If the system contains an isocratic pump go on to Step 16.

Note: Microbore (MB) and Standard Bore (SB) pump configurations use different sequences for the Gradient Accuracy Test. GP40 and GP50 microbore pumps should be qualified using the GP40_GP50_MB_Grad_Accuracy sequences. GS50 and standard bore GP50 pumps should be qualified using the GP50_GS50_SB_Grad_Accuracy sequences. Standard bore GP40 pumps should be qualified using the GP40_SB_Grad_Accuracy sequences.

16. If the system to be qualified contains an ED40 or ED50 Detector used as an Amperometry Detector, select the **Int_Amp_Noise_Drift** and/or the **DC_Amp_Noise_Drift** sequence template. Copy the template into the **Validation** directory. These tests should only be performed if the ED40 or ED50 Detector is used in the Integrated Amperometry or DC Amperometry mode.
17. If the system to be qualified contains and EG40 Eluent Generator, select the **EG40 Current Qual** sequence template and copy it into the **Validation** directory.
18. Lastly, copy the report file into the **Validation** directory. Select the **OQ_UV.rdf** or the **PQ_UV.rdf** for qualification with an AD20, AD25 or PDA-100 detector. Select the **OQ_ECD.rdf** or the **PQ_ECD.rdf** file for qualification with a Conductivity Detector.
19. Open the **Validation** directory. Right click on the first sequence template and select **Properties**. Select the timebase that will be qualified in the **Timebase** entry field. In the **Preferred Report and Channel** section, select the report (rdf) file newly copied into the directory of the timebase. Confirm that the correct detector channel is selected.
20. Repeat step 19 for all other OQ/PQ sequence templates in the **Validation** directory.

2.7 OQ/PQ Report Set Up

All relevant information such as the serial number, module type, moduleware version and customer information must be entered into the SPECIFICATION sheet of the OQ and PQ report file.

1. Open the OQ.rdf or PQ.rdf file copied into the validation directory in Section 2.6. Locate the row of sheet names at the bottom of the report. Click on the first sheet of the report titled **SPECIFICATION**. Go to **Edit** and turn **Layout Mode** on.
2. In the **Instruments** section of the **SPECIFICATION** sheet, enter the **Model**, **Supplier's Name** (if other than Dionex), **Serial Number** and **Moduleware Version** for the pump, detector, autosampler, software and LC module configured on the IC system. Type N/A in the boxes for any modules in the SPECIFICATION sheet that were not used.
3. In the Additional Information section located on page 2 of the **SPECIFICATION** sheet, enter the **Customer Name/Company** and **Qualification Executor/Company** information.
4. Right click anywhere on the sheet and select **Save Report Definition**. Click on **Edit** and turn **Layout Mode** off. Close the report file.

3 Performing the Qualification

3.1 When to Perform Operational Qualification

Operational Qualification (OQ) should be performed following the initial installation of an IC system and on any new individual module installed on an existing system.

3.2 When to Perform Performance Qualification

Performance Qualification (PQ) can be performed at any time after the initial installation and OQ. Dionex recommends that a performance qualification be performed on IC systems every 6 months. PQ should be performed only after preventive maintenance has been completed on the system.

3.3 OQ/PQ Test Order

The OQ/PQ qualification tests are performed in the following order:

- A. Wavelength Accuracy Test for Absorbance Detectors
- B. Detector Noise and Drift Test
- C. Injector Precision Test
- D. Detector Linearity Test
- E. Injector Linearity Test - *AS50 only*
- F. Carry Over Test
- G. Sample Prep Precision Test - *AS50 with Sample Prep Option only*
- H. Gradient Accuracy Test*
- I. DC Amperometry Baseline Noise and Drift*
- J. Integrated Amperometry Baseline Noise and Drift*
- K. EG40 Current Qualification Test*

* Perform this test only if module is configured and/or used in this mode on the timebase to be qualified.

3.4 What to Do if a Test Fails

If at any point the instrument results do not pass within Dionex specified limits, the instrument must be serviced. See Section 5 for troubleshooting assistance. When the repair is complete, rerun all tests from the original point of failure on in the sequence. For example, if the instrument passes Tests A, B, and C and fails Test D, make the appropriate repairs and then repeat tests D - G only.

3.5 Starting the Sequences

The OQ and PQ sequences are system-specific test procedures. If different module options or configurations are used, the program files located in the standard sequences must be modified accordingly.

3.5.1 Wavelength Accuracy

3.5.1.1 AD25 and PDA - 100

The flow cell must be free from colored reagents before this test is run. If there is colored reagent in the flow cell (for example, PAR used in transition metal analysis), the flow cell must be thoroughly flushed prior to running the test.

1. Select 100% water and a flow rate of 1.00 mL/min. Turn the pump on and allow the detector to equilibrate.

3.5.1.1.1 Wavelength Accuracy from a Wellness Panel

1. If PeakNet 6.11 is installed, open the **Dionex_pump_absorbance_wellness** panel from the \Panels\Wellness directory to perform the wavelength accuracy test. This panel will work for the AD25 and PDA -100 detectors.
2. **Note:** The AD25 detector must be configured for System Wellness to use the Wellness Panel. If Wellness is not configured, perform the Wavelength accuracy test from the AD25 front panel using the steps in Section 3.5.1.1.2.
3. Go to the Diagnostic Tests section of the panel, then to the Wavelength Verification section and select the **verify** button.
4. After approximately two minutes, the Wellness panel will display the theoretical and measured wavelengths, and the delta between the theoretical and measured values.
5. Open the **OQ_UV.rdf** or **PQ_UV.rdf** report file, depending on the qualification being performed. Locate the row of sheet names at the bottom of the report. Click on the **DET_AD25_PDA_WAVE** sheet.
6. Go to **Edit** and turn **Layout Mode** on. Enter today's date in the **Additional Information** section. Enter the **Measured** wavelength values into the table on the sheet. The deviation and pass/fail results will be calculated automatically. Print the sheet. Click on **Edit** and turn **Layout Mode** off. Close the report file.
7. If any part of the test does not pass, refer to Section 5 for troubleshooting assistance.

3.5.1.1.2 Wavelength Accuracy from the AD25 Front Panel

This test can be performed on the AD25 only.

1. From the AD25 front panel, access the Main Menu and press (8) to enter the Diagnostic Menu.
2. From the AD25 Diagnostic Menu, press (5) for Wavelength Verify.
3. Select **Run** and press **Enter**.

4. After approximately two minutes, the front panel will display the theoretical and measured wavelengths, and the delta between the theoretical and measured values.
5. Open the **OQ_UV.rdf** or **PQ_UV.rdf** report file, depending on the qualification being performed. Locate the row of sheet names at the bottom of the report. Click on the **DET_AD25_PDA_WAVE** sheet.
6. Go to **Edit** and turn **Layout Mode** on. Enter today's date in the **Additional Information** section. Enter the **Measured** wavelength values into the table on the sheet. The deviation and pass/fail results will be calculated automatically. Print the sheet. Click on **Edit** and turn **Layout Mode** off. Close the report file.
7. If any part of the test does not pass, refer to Section 5 for troubleshooting assistance.

3.5.1.2 AD20

1. Remove the AD20 Sample cell from the photodiode and disconnect the beamsplitter cable.
2. Install the AD20 Holmium Oxide Test cell (Dionex P/N 046172) onto the photodiode and connect the beamsplitter cable.
3. Make sure the UV and Vis lamps are set to HI. Set the detector wavelength to 354 nm and then press **Enter** on the front panel. Allow the detector to equilibrate with the test cell installed for at least 10 minutes before performing the Wavelength Accuracy test.
4. Open the **OQ_UV.rdf** or **PQ_UV.rdf** report file, depending on the qualification being performed. Locate the row of sheet names at the bottom of the report. Click on the **DET_AD20_WAVE** sheet. Print the sheet. This will allow for manual recording of test results. Test can also be directly entered into the report (see Step 10). Close the report.
5. From the front panel of the AD20 detector, press the **Offset** button to zero the detector reading. Allow the detector to equilibrate for approximately 20 seconds. Record the offset absorbance value as the measured absorbance for 354 nm. This value should be approximately 0.0000 AU.
6. Manually change the wavelength setting on the front panel of the AD20 to 355 nm and press **Enter** on the front panel. Allow the detector to equilibrate for approximately 20 seconds. Record the resulting absorbance value for 355 nm.
7. Repeat Step 6 for the remaining 13 wavelength settings in the **DET_AD20_WAVE** sheet.
8. After all wavelength absorbance values have been recorded, find the wavelength setting that produced the maximum absorbance value.
9. Record the wavelength setting that produced the maximum absorbance value.
10. Open the **OQ_UV.rdf** or **PQ_UV.rdf** report file and go to the **DET_AD20_WAVE** sheet. Go to **Edit** and turn **Layout Mode** on. Enter today's date in the **Additional Information** section. Enter all previously recorded values in the table on the sheet. The deviation and pass/fail results will be calculated automatically. Print the sheet. Click on **Edit** and turn **Layout Mode** off. Close the report file.
11. If any part of the test does not pass, refer to Section 5 for troubleshooting assistance.

3.5.2 Performing the Main Qualification

The main part of the OQ/PQ consists of the tests in the Detector_Pump sequence. They include:

Detector Noise and Drift (Absorbance and Conductivity detectors)

Injector Precision

Detector Linearity

Injector Linearity (AS50 sequences only)

Carry Over

Sample Prep Precision (AS50 SP sequences only)

This part of the qualification requires that an AD20, AD25, PDA-100 or Conductivity Detector be configured on the system being qualified. If an Amperometry Detector is the only detector configured on the system, go on to Section 3.5.4 for DC Amperometry qualification or Section 3.5.5 for Integrated Amperometry qualification.

If the qualification will be performed using a manual injection valve, use full loop injections and overfill the loop a minimum of 3 times. Use the AS40 sequence templates. Inject standards at the beginning of each sample in the sequence. There will be a 2.5 minute delay before the sample is injected.

Note: AD20 and AD25 sequences contain both 254 nm and 520 nm Noise and Drift Tests. If you do not want to perform the Noise and Drift tests at both wavelengths, change the sequence status line of the test you do not want to run to **Finished**, then save the sequence. The PDA - 100 sequences collect data at both wavelengths simultaneously.

1. The IC must be set up and running as described in Section 2.5.
2. From the Browser, go to the Panel directory and open the panel that best matches your system configuration. Go to **Control** and select the timebase to be qualified.
3. Go to **Batch** and **Edit**. **Remove** any existing sequences. Click on **Add** and select the Detector_Pump sequence chosen in Section 2.6.
4. Click on **Ready Check** to confirm the program and method files contained in the sequence are compatible with the configured system. Make modifications to the program files if necessary.

Note: If an EG40 is configured on the system, the following message will occur at the Ready Check: **"Flow rate not specified at time 0.0, using a flow setting 1.00."** This refers to the EG40 flow rate, not the pump flow rate. It is not necessary to take any action since the EG40 is not used for this part of the qualification. Proceed with starting the batch.

5. **Start** the Batch. Run times are as follows:

AS40 Detector_Pump sequences have a run time of approximately 1 hr and 10 mins.

AS50 CC/TC Detector_Pump sequences have a run time of approximately 2 hrs, 10 mins.

AS50 CC/TC SP Detector_Pump sequences have a run time of approximately 2 hrs, 40 mins.

At the end of the sequence, the system will automatically stop.

6. When the Batch is finished, go to the Browser and right click on the Detector_Pump sequence. Go to **Batch Report**. Confirm that the correct report (rdf) file is selected and **Print** the report.
7. Review the report for Pass/Fail results. If any test does not pass, make the appropriate instrument repair/adjustments and rerun all tests from the original point of failure on. Refer to Section 5 for troubleshooting assistance.

3.5.3 Gradient Accuracy

The Gradient Accuracy test only needs to be performed on gradient pumps. If the gradient pump is used only in isocratic mode, this test does not need to be run.

The Gradient Accuracy test can be run on standard bore and microbore GP40 and GP50 pumps, and on all GS50 pumps. GP40 and GP50 microbore pumps should be qualified using the GP40_GP50_MB_Grad_Accuracy sequences. GS50 and standard bore GP50 pumps should be qualified using the GP50_GS50_SB_Grad_Accuracy sequences. Standard bore GP40 pumps should be qualified using the GP40_SB_Grad_Accuracy sequences.

1. Make sure all pump lines are primed well, according to the system preparation in Section 2.5.
2. Turn the pump flow off.
3. Remove the Backpressure Coil and replace it with approximately 30" of yellow PEEK tubing.
4. Set the pump flow rate to 1.00 mL minute. Monitor the system backpressure and confirm the reading is between 1500 and 2000 psi. If the backpressure reading falls outside of this range, adjust the length of the yellow PEEK tubing accordingly.
5. Allow the pump to equilibrate for approximately 10 minutes.
6. From the Browser, go to the Panel directory and open the panel that best matches your system configuration. Go to **Control** and select the timebase to be qualified.
7. Go to **Batch** and **Edit**. Click on **Add** and select the Grad Accuracy sequence chosen in Section 2.6. The sequence templates are different for GP40 and GP50/GS50 pumps.
8. Click on **Ready Check** to confirm the program and method files contained in the sequence are compatible with the configured system. Make modifications to the program files if necessary.
9. **Start** the Batch. Run times are as follows:

GP40 Grad Accuracy sequences have a run time of approximately 37 mins.
GP50_GS50 sequences have a run time of approximately 62 mins.
10. At the end of the sequence, the system will automatically stop.

11. When the Batch is finished, go to the Browser and right click on the Grad Accuracy sequence. Go to **Batch Report**. Confirm that the correct report (rdf) file is selected and **Print** the report.
12. Review the report for Pass/Fail results. If any part of the test does not pass, make the appropriate repair/adjustments and rerun the test. Refer to Section 5 for troubleshooting assistance.

3.5.4 DC Amperometry Noise and Drift

This test should only be performed if the ED40 or ED50 Detector is used in the DC Amperometry mode.

1. Disconnect the Amperometry cell cable from the ED40/ED50 Detector.
2. Plug the ED/CD Validation Test Cell Amperometry cable into the ED40/ED50 Amperometry Cell connector.
3. From the front panel of the detector, turn the cell on.
4. Let the detector equilibrate with the validation cell for approximately five minutes.
5. From the Browser, go to the Panel directory and select the **ED40_DC_Amp panel**. This panel will work with either the ED40 or ED50 detector. Go to **Control** and select the timebase configured with the detector to be qualified.
6. Go to **Batch** and **Edit**. Click on **Add** and select the **DC_Amp_Noise_Drift** sequence chosen in Section 2.6.
7. Click on **Ready Check** to confirm the program and method files contained in the sequence are compatible with the configured system. Make modifications to the program files if necessary.
8. **Start** the Batch. The run time is 20 minutes. The pump flow for this sequence is 0.00 mL/min.
9. When the Batch is finished, go to the Browser and right click on the **DC_Amp_Noise_Drift** sequence. Go to **Batch Report**. Confirm that the OQ_ECD or PQ_ECD report (rdf) file is selected and **Print** the report.
10. Review the report for Pass/Fail results. If any part of the test does not pass, make the appropriate repair/adjustments and rerun the test. Refer to Section 5 for troubleshooting assistance.

3.5.5 Integrated Amperometry Noise and Drift

This test should only be performed if the ED40 or ED50 Detector is used in the Integrated Amperometry mode.

1. Disconnect the Amperometry cell cable from the ED40/ED50 Detector.
2. Plug the ED/CD Validation Test Cell Amperometry cable into the ED40/ED50 Amperometry Cell connector.
3. From the front panel of the detector, turn the cell on.
4. Let the detector equilibrate with the validation cell for approximately five minutes.
5. From the Browser, go to the Panel directory and select the **ED40_Integrated_Amp panel**. This panel will work with either the ED40 or ED50 detector. Go to **Control** and select the timebase configured with the detector to be qualified.
6. Go to **Batch** and **Edit**. Click on **Add** and select the **Int_Amp_Noise_Drift** sequence chosen in Section 2.6.
7. Click on **Ready Check** to confirm the program and method files contained in the sequence are compatible with the configured system. Make modifications to the program files if necessary.
8. **Start** the Batch. The run time is 20 minutes. The pump flow for this sequence is 0.00 mL/min.
9. When the Batch is finished, go to the Browser and right click on the **Int_Amp_Noise_Drift** sequence. Go to **Batch Report**. Confirm that the OQ_ECD or PQ_ECD report (rdf) file is selected and **Print** the report.
10. Review the report for Pass/Fail results. If any part of the test does not pass, make the appropriate repair/adjustments and rerun the test. Refer to Section 5 for troubleshooting assistance.

3.5.6 EG40 Current Qualification

The EG40 Current Qualification test can be used for both KOH and MSA EG40 Cartridges.

1. Disconnect the Eluent Generator Cartridge cable from the EG40.
2. Plug the EG40 Validation Test Cell cable into the EG40 Cartridge connector.

Note: This test requires EG40 Validation Test Cell Rev 02 or higher. To modify a Rev 01 cell to use for the current test, see Section 8.

3. Turn on the Voltmeter and set the display to mA. This test requires mA readings from 0.00 to 220 mA.

4. Connect the positive lead of the voltmeter to the Test Cell red wire connection closest to the EG40 cartridge connection.
5. Connect negative lead of the voltmeter to the other Test Cell red wire connection.
6. From the Browser, go to the Panel directory and select the pump_EG40_detector panel appropriate for the system being qualified. Go to **Control** and select the timebase configured with the EG40 to be qualified.
7. From the EG40 section of the panel, select the cartridge type to be used. The panel will recognize previously configured KOH and MSA cartridges even though the EG40 Validation cell is plugged into the cartridge connector.
8. Go to **Batch** and **Edit**. Click on **Add** and select the **EG40 Current Qual** sequence chosen in Section 2.6.
9. Click on **Ready Check** to confirm the program and method files contained in the sequence are compatible with the configured system. Make modifications to the program files if necessary. The pump flow for this sequence is set to 1.00 mL/min to keep the pressure above the minimum pressure limit required for the EG40.
10. **Start** the Batch. This sequence is an EG40 step gradient. Each step is held for one minute. There is a delay of 0.4 seconds between when the EG40 program steps to the next concentration and when the voltmeter registers the new milliamp reading. Record the milliamp reading for each step. These values will be entered into the **EG40_CUR_KOH** or **EG40_CUR_MSA** sheets after the run is complete. The run time is 16 minutes. Data acquisition is off during this run.
11. When the run is complete, open the **OQ_UV.rdf** or **PQ_UV.rdf** report file, depending on the qualification being performed. Locate the row of sheet names at the bottom of the report. Click on the **EG40_CUR_KOH** or **EG40_CUR_MSA** sheet, depending on the type of cartridge to be qualified.
12. Go to **Edit** and turn **Layout Mode** on. Enter today's date in the **Additional Information** section. Enter the milliamp reading for each step in the **EG40_CUR_KOH** or **EG40_CUR_MSA** sheet. The deviation and pass/fail results will be calculated automatically. Print the sheet, click on **Edit** and turn **Layout Mode** off. Close the report.
13. If any part of the test does not pass, make the appropriate repair/adjustments and rerun the test. Refer to Section 5 for troubleshooting assistance.
14. To qualify another cartridge, select the new cartridge type in the control panel and repeat steps 7 – 12.

3.6 OQ/PQ Completion

3.6.1 Customer Review

1. The results should be reviewed by the instrument owner/user. If the qualification is accepted, both the Customer and Qualification Executor need to sign the printed reports.
2. Fill out the Validation Stickers and Certificates of Functional Performance for the qualified modules.
3. Place the reports and certificates into the binder provided in the OQ/PQ kit and leave it with the customer.

3.6.2 Returning the System to the Original Configuration

1. If the Signal factor was modified in the system preparation, it must be set back to its original value. Open **Server Configuration**. Double click on the detector/s that were qualified to view the properties page. Select the **Signals** page and change the factor back to the original value. Exit the properties page and save the configuration.
2. Disconnect the Backpressure Coil and/or yellow PEEK tubing and restore the system plumbing to its original configuration with the appropriate columns, suppressors, and EG40 in line.
3. Install the original sample loop back onto the injection valve. For systems using an AS50 Autosampler, go to the **Main Menu** on the front panel of the autosampler. Access the **Module Set Up** screen and then go to **Plumbing Configuration**. Change the sample loop size back to its original value.
4. Reconnect the original eluents to the pump and prime each line individually.
5. Flush the system with an eluent recommended by the customer and confirm that there are no leaks.

4 Explanation of Test Procedures

4.1 Wavelength Accuracy of Absorbance Detectors

4.1.1 Theory

Wavelength accuracy is evaluated using a holmium oxide filter. Three holmium oxide maxima are compared with the theoretical maxima at 360.9 nm, 418.7 nm and 536.2 nm for AD25 and PDA – 100 detectors. Absorbance readings are measured with a holmium oxide filter at varying wavelength settings for the AD20 wavelength accuracy test.

4.1.2 Performing and Evaluating the Test

The wavelength accuracy test is triggered from the front panel of the AD25 or through a PeakNet™ Wellness panel for either the AD25 or PDA-100. When the verification is complete, test results are displayed either on the front panel of the AD25 or in the Wellness panel. Measured wavelengths, deltas and pass/fail results are recorded manually. The test is performed from the front panel of the AD20 by setting wavelengths manually and recording the resulting absorbance values.

4.2 Baseline Noise and Drift

4.2.1 Theory

A 20-minute detector baseline is collected using water as the eluent at a flow rate of 1.00 mL/min. For absorbance detectors, data is collected at 254 nm, 520 nm, or both.

4.2.2 Performing the Test and Evaluating the Results

Noise and drift evaluation for all IC detectors is the first test in the Detector_Pump sequence and is called Noise and Drift. This is Test B in the OQ/PQ test order.

Noise is calculated by splitting the 20-minute baseline signal into 20 one-minute intervals. The noise is the distance between two parallel lines through the measured minimum and maximum values and the regression line. The measured noise values are then averaged to establish the final value. To calculate drift, PeakNet™ uses a method of least squares and calculates a regression line from the measured values. The slope of the curve indicates the drift of the measured signal, and the amount of the slope indicates the amount of drift.

4.3 Injector Precision

4.3.1 Theory

Injector precision is determined by making six replicate injections of a standard using a fixed sample volume. A 40 ppm caffeine standard is injected when using absorbance detection at a wavelength of 272 nm. A 50 ppm nitrate standard is injected when using conductivity detection. Water is used as the eluent at a flow rate of 1.00 mL/min.

4.3.2 Performing and Evaluating the Test

Injector precision evaluation is the second set of samples in the Detector_Pump sequence and the samples are called Inj Precision_1 through Inj Precision_6. This is Test C in the OQ/PQ test order. The concentration of the nitrate or caffeine standard and the sample volume injected is included in the QNT file used for injector linearity.

Precision is determined by calculating the relative standard deviation of the peak areas for the six injections.

4.4 Detector Linearity

4.4.1 Theory

Detector linearity for absorbance detectors is determined at 272 nm using a sequence of five caffeine standard injections at 15, 40, 60, 80 and 100 ppm. Linearity for conductivity detectors uses a sequence of five nitrate standard injections at 5, 10, 25, 50, and 100 ppm. Water is used as the eluent at a flow rate of 1.00 mL/min.

4.4.2 Performing the Test and Evaluating the Results

Detector linearity evaluation is the third test in the Detector_Pump sequence and the samples are called Detector linearity_1 through Detector linearity_5. This is Test D in the OQ/PQ test order. The concentrations for the nitrate and caffeine standards are included in the QNT file used for detector linearity.

PeakNet™ calculates the peak height for each standard. A plot of peak height vs. concentration is generated and the regression line is determined. The correlation coefficient and the relative standard deviation of this line indicate linearity.

4.5 Injection Volume Linearity

4.5.1 Theory

Injection volume linearity is determined by making five injections of a standard using injection volumes of 10, 20, 30, 40 and 50 μL . A 60 ppm caffeine standard is injected when using absorbance detection at a wavelength of 272 nm. A 100 ppm nitrate standard is injected when using a conductivity detection. Water is used as the eluent at a flow rate of 1.00 mL/min. This test is only performed on autosamplers with variable injection volume capability.

4.5.2 Performing the Test and Evaluating the Results

Injection volume linearity evaluation is the fourth test in the Detector_Pump sequence and the samples are called Injector linearity_1 through Injector linearity_5. This is Test E in the OQ/PQ test order. The concentrations for the nitrate and caffeine standard and the injection volumes are included in the QNT file used for detector linearity.

PeakNet™ calculates the peak area for each injection. A plot of peak area vs. injection volume is generated and the regression line is determined. The correlation coefficient and the relative standard deviation of this line indicate linearity.

4.6 Sample Carry-over

4.6.1 Theory

Sample carry-over is determined by injecting a highly concentrated sample followed by a blank injection. The blank in this case is water. Nitrate at a concentration of 1000 ppm is injected when this test is performed with a conductivity detector. Caffeine at a concentration of 1000 ppm is injected when this test is performed with an absorbance detector. Ideally, the chromatogram of the blank injection should not contain any peaks. If, however, a peak is detected, the ratio of the peak areas is compared. Water is used as the eluent at a flow rate of 1.00 mL/min.

4.6.2 Performing the Test and Evaluating the Results

Sample carry-over evaluation is the fifth test in the Detector_Pump sequence and the samples are called Carry over and Water. This is Test F in the OQ/PQ test order.

A peak area ratio is calculated for the two injections determine sample carry-over.

4.7 Sample Prep Precision

4.7.1 Theory

Sample Prep precision is determined by making three separate injections using a fixed injection volume. Each of the three standards is prepared individually using the autosampler's sample preparation option. The autosampler prepares the three samples by first transferring 15 μL of a 1000 ppm caffeine or nitrate standard into three empty sample vials and then adding 1485 μL of water to each. Injections are then made from each vial. Water is used as the eluent at a flow rate of 1.00 mL/min. If an absorbance detector is used the wavelength is set to 272 nm. This test is only performed on autosamplers with the sample prep option.

4.7.2 Performing the Test and Evaluating the Results

Sample prep precision evaluation is the last set of samples in the Detector_Pump sequence and the samples are called Sample Prep_1 through Sample Prep_3. This is Test G in the OQ/PQ test order.

Precision is determined by calculating the relative standard deviation of the peak areas for the three injections.

4.8 Gradient Accuracy

4.8.1 Theory

Gradient accuracy is determined by performing a step gradient using varying compositions of water and either 0.3% Acetone or 1 mM KCL. When performing the test with an absorbance detector, water is used for eluents A and C and 0.3% Acetone is used for eluents B and D with a wavelength setting of 254 nm. When performing the test with a conductivity detector, water is used for eluents A and C and 1 mM KCl is used for eluents B and D. Gradient composition is measured at 0, 20, 50, 80, and 100% of eluents B and D at a flow rate of 1.00 mL/min with a backpressure of 1500 -2000 psi. Step noise (ripple) is measured for each step in the gradient. No injection is required.

4.8.2 Performing the Test and Evaluating the Results

Gradient Accuracy evaluation is contained in the GP40_Grad Accuracy or GP50_GS50 Grad Accuracy sequences. This is Test H in the OQ/PQ test order.

The measured step height is compared to the theoretical step height and a percent deviation is calculated. Step noise is evaluated for all steps. Each step has a defined one-minute interval where noise is evaluated. Based on a method of least squares, PeakNet™ calculates a regression line from the measured values of each interval. The noise is the distance between two parallel lines through the measured minimum and maximum values and the regression line.

4.9 DC Amperometry Noise and Drift

4.9.1 Theory

A 20-minute detector baseline is collected using an amperometry validation test cell.

4.9.2 Performing the Test and Evaluating the Results

Noise and drift evaluation for the DC Amperometry detector is contained in the DC_Amp_Noise_Drift sequence. This is Test I in the OQ/PQ test order.

Noise is calculated by splitting the 20-minute baseline signal into 20 one-minute intervals. The noise is the distance between two parallel lines through the measured minimum and maximum values and the regression line. The measured noise values are then averaged to establish the final value. To calculate drift, PeakNet™ uses a method of least squares and calculates a regression line from the measured values. The slope of the curve indicates the drift of the measured signal, and the amount of the slope indicates the amount of drift.

4.10 Integrated Amperometry Noise and Drift

4.10.1 Theory

A 20-minute detector baseline is collected using an amperometry validation test cell.

4.10.2 Performing the Test and Evaluating the Results

Noise and drift evaluation for the Integrated Amperometry detector is contained in the Int_Amp_Noise_Drift sequence. This is Test J in the OQ/PQ test order.

Noise is calculated by splitting the 20-minute baseline signal into 20 one-minute intervals. The noise is the distance between two parallel lines through the measured minimum and maximum values and the regression line. The measured noise values are then averaged to establish the final value. To calculate drift, PeakNet™ uses a method of least squares and calculates a regression line from the measured values. The slope of the curve indicates the drift of the measured signal, and the amount of the slope indicates the amount of drift.

4.11 EG40 Current Qualification

4.11.1 Theory

EG40 Current Qualification is determined by performing an EG40 step gradient using varying eluent concentration set points. The EG40 applies a specific current to the EluGen cartridge depending on the concentration of eluent and type of cartridge selected. An EG40 validation test cell is used in place of the EluGen cartridge and the applied current for each step is measured using a voltmeter that is connected to the test cell.

4.11.2 Performing the Test and Evaluating the Results

The EG40 current qualification test is contained in the EG40_Current_Qual sequence. This is Test K in the OQ/PQ test order. Current readings in mA are displayed on the voltmeter for each concentration step in the gradient. Values for each step and pass/fail results are manually recorded.

5 Troubleshooting

5.1 General Hints

If the system pressure is above 1500 psi after connecting the PEEK Backpressure coil with a flow rate of 1.00 mL/min, the coil has a blockage. Replace the coil.

Dips before or after the main nitrate or caffeine peak indicate extra dead volume in the system. This is caused by small voids between the end of a liquid line and the bottom of the port it is installed in. One by one disconnect each line from the injection valve to the detector cell. Push the ferrule and fitting up and away from the end of the tubing. Make sure the tubing is cut straight. If it is not, recut the end. When reinstalling all liquid connections, first insert the tubing all the way into the port until it touches the bottom. Hold the tubing in this position while inserting and tightening ferrule and fitting down into the port.

5.2 Failure of Individual Tests

5.2.1 Absorbance Detector

Test	Reason	Action
Wavelength Accuracy	Optics not calibrated	Perform wavelength calibration and then repeat the wavelength accuracy test
	Air bubbles in flow cell	Flush out flow cell
	Colored reagent in flow cell	Flush out flow cell
Baseline Noise	Poor water quality	Replace with an alternate water source
	Lamp is reaching maximum lifetime	Replace lamp
	Air bubbles in flow cell	Flush out flow cell
Drift	Detector not warmed up	Allow recommended warm up period (refer to test procedure)
	System not equilibrated	Allow longer flush/equilibration time
	Lamp defective	Replace lamp
Lamp Intensity (PDA-100 only)	Lamp is reaching maximum lifetime	Replace lamp
	Flow cell inserted incorrectly	Reinstall flow cell
Detector Linearity	Lamp is reaching maximum lifetime	Replace lamp
	Concentration of standards is incorrect	Remake standards
	Liquid leak	Find source of leak

5.2.2 Conductivity Detector Tests

Test	Reason	Action
Baseline Noise	Poor water quality Air bubbles in flow cell Cell temperature regulation failure	Replace with an alternate water source Flush out flow cell Replace cell
Drift	System not equilibrated Lab temperature not stable	Allow longer flush/equilibration time Regulate temperature (block vents, close windows, etc.)
Detector Linearity	Concentration of standards is incorrect Liquid leak Injection volume not consistent (manual injections only)	Remake standards Find source of leak Overfill sample loop at least 3x when loading.

5.2.3 Autosampler

Test	Reason	Action
Injector Precision	Air injected from vial or from sample syringe when doing manual injections	Too little sample volume in vial - refill vial or flush bubbles out of injection syringe.
	Air bubbles in AS50 sample syringe	Flush out syringe
	Liquid leak in autosampler or injection valve	Find source of leak
Injection Volume Linearity	Detector failure	See Detector linearity troubleshooting section
	Defective sample syringe	Replace sample syringe
Carry Over	Contaminated injection valve	Replace/clean injection valve rotor and stator face
	Contamination in autosampler	Flush autosampler liquid lines completely
	Contamination in sample vial or in syringe used for manual injections	Use new sample vial for water blank, rinse vial thoroughly before using, flush syringe or use a new one
Sample Prep Precision	Air injected from vial Liquid leak	Too little sample volume in vial, refill vial.
	Air bubbles in AS50 sample syringe or prep syringe	Flush out syringe
	Defective prep syringe	Replace prep syringe
	Liquid leak in autosampler or injection valve	Find source of leak

5.2.4 Pump

Test	Reason	Action
Gradient Step Accuracy	Air in system	Prime pump - prime each line individually for 3 - 5 minutes, flush system
	System not equilibrated	Allow longer flush/equilibration time.
	Old eluent in vacuum degas chamber	Prime pump - prime each line individually for 3 - 5 minutes
	Incorrect eluent concentration for eluents B and D	Remake 0.3% Acetone or 1mM KCL solutions
	GP40 pump eluents not pressurized	Pressurize eluent bottles
Step Noise (ripple)	GP50/GS50 eluent bottles not vented	Vent eluent bottles
	Air in system	Prime pump, flush system
Step Noise (ripple)	Air in flow cell	Flush flow cell
	System not equilibrated	Allow longer flush/equilibration time
	Check valve or piston seal failure	Replace check valve/piston seal

5.2.5 DC and Integrated Amperometry Detector Tests

Test	Reason	Action
Baseline Noise	Validation cell not connected properly	Reconnect test cell cable to amperometry cell connection on SP board
	Validation cell defective	Replace validation test cell
Drift	Detector not equilibrated	Allow longer equilibration time
	Lab temperature not stable	Regulate temperature (block vents, close windows, etc.)

5.2.6 EG40

Test	Reason	Action
Current	Validation cell not connected properly	Reconnect test cell cable to the EG40, confirm voltmeter leads are making contact with test cell wires
	Voltmeter on wrong setting	Voltmeter should be set to read current (mA). Range should be set to view 0.00 to 210.00 mA
	Validation cell defective	Replace validation test cell

6 Modifying an EG40 Validation Test Cell

If using an EG40 Rev 01 Validation Test Cell, the cable needs to be modified to perform the EG40 Current Qualification Test with PeakNet 6. Rev 02 Validation Test Cells have been modified for use with this test already.

1. Using a razor blade carefully cut the heat shrink tubing on the test cell cable lengthwise. Make a cut approximately 2 inches long in the cable mid-section. This will expose a black and red wire.
2. Cut the red wire and strip off about a ¼ inch of the red insulation, exposing the bare wire. Connect each end of the bare wire to the voltmeter as directed in Section 3.5.6.

7 Example Reports

Following are two example PQ Reports. The PQ_UV report is for a system configured with an AD25, GS50, and an AS50 with a CC. An example AD20 Wavelength Accuracy report is also included at the end of the PQ_UV report.

The PQ_ECD report is for a system configured with an ED40 (conductivity and integrated amperometry modes), GP40, EG40 and an AS50 with a CC. The OQ and PQ reports are the same except the test limits differ according to the table in Section 1.3.



Performance Qualification

Instruments

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GS50	Dionex	00071255	1.00
LC Module or LC Oven	CC	Dionex	99125564	N/A
Absorbance Detector	AD25	Dionex	00012233	1.02
Chromleon Datasystem	6.20 Build 527 (Demo-Installation)	Dionex	12345678	N/A
Autosampler	AS50	Dionex	00052145	1.07

Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
Backpressure Tubing	0.003" ID PEEK, 50"
UV Sample 1	Caffeine, 15 ppm
UV Sample 2	Caffeine, 40 ppm
UV Sample 3	Caffeine, 60 ppm
UV Sample 4	Caffeine, 80 ppm
UV Sample 5	Caffeine, 100 ppm
UV Sample 6	Caffeine, 1000 ppm
Blank	Water
Eluent A	Water
Autosampler Reservoir A	Water
AD20 Holumium Oxide Test Cell	

Customer Signature

Qualification Executor

Date

Customer signature indicates that all information in the following reports has been reviewed and accepted.

Limits

Test	Customized Limits	Dionex Recommended Limits
Noise - AD20 UV and Vis (µAU)	50.000	50.000
Drift - AD20 UV and Vis (µAU/hr)	200.000	200.000
Noise - AD25 UV and Vis (µAU)	40.000	40.000
Drift - AD25 UV and Vis (µAU/hr)	200.000	200.000
Noise - PDA - 100 UV and Vis (mAU)	0.100	0.100
Drift - PDA - 100 UV and Vis (mAU/hr)	1.000	1.000
AD25, PDA - 100 Wavelength Accuracy (nm, +/-)	2	2
AD20 Wavelength Accuracy (nm, +/-)	3	3
Detector Linearity (Corr.)	0.999	0.999
Detector Linearity (%RSD)	5.000	5.000
Injector Precision (Area %RSD)	1.000	1.000
Carry Over (Area %)	0.100	0.100
Injector Linearity (Corr.)	0.999	0.999
Injector Linearity (%RSD)	1.000	1.000
Sample Prep Precision (%RSD)	2.000	2.000
Standard Bore Gradient Pump Step Accuracy %	2.000	2.000
Standard Bore Pump Noise	0.20%	0.20%
Microbore Gradient Pump Step Accuracy %	2.000	2.000
Microbore Pump Noise	0.50%	0.50%

Additional Information

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
 Qualification Executor/Company: **Mr. Smith / Dionex** Next Qualification: **Feb-01**

Customer Signature _____ Qualification Executor _____ Date _____



Performance Qualification

Wavelength Accuracy

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Absorbance Detector	AD25	Dionex	00012233	1.02
Chromeleon Datasystem	6.20 Build 527 (Demo-Installation)	Dionex	12345678	N/A

• Accessories

<i>Name</i>	<i>Description</i>
None	

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company**

Date: **08-28-00**

Qualification

Executor/Company: **Mr. Smith / Dionex**

Next Qualification: **Feb-01**

• Data for the Wavelength Accuracy Test

<i>Theoretical Wavelength (nm)</i>	<i>Measured Wavelength (nm)</i>	<i>Delta (nm)</i>	<i>+/- Limit (nm)</i>	<i>Result</i>
360.9	361.0	0.1	2	PASS
418.7	418.9	0.2	2	PASS
536.2	536.3	0.1	2	PASS

Customer Signature

Qualification Executor

Date



Performance Qualification

AD25 Detector Noise and Drift: 520 nm

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GS50	Dionex	00071255	1.00
Absorbance Detector	AD25	Dionex	00012233	1.02

• Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
Eluent	Water

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
Qualification Executor/Company: **Mr. Smith / Dionex** Next Qualification Due: **Feb-01**

• Test Results Summary

<i>Test</i>	<i>Result</i>
Noise	PASS
Drift	PASS

Customer Signature

Qualification Executor

Date

• **Data for Detector Noise**

Segment No.	Noise, μAU
1	16.294
2	19.354
3	9.218
4	13.404
5	13.488
6	11.191
7	15.872
8	17.090
9	15.184
10	14.510
11	11.528
12	12.572
13	14.657
14	11.976
15	12.794
16	11.333
17	11.973
18	13.820
19	15.897
20	11.623
Average (μAU)	13.958
Limit (μAu)	40.000
Result	PASS

• **Data for Detector Drift**


20 Minute Drift, μAU	Drift, μAu/hr	Limit, μAu/hr	Result
0.047	0.140	150.000	PASS

Customer Signature

Qualification Executor

Date

• **Chromatogram of Detector Noise and Drift**



Customer Signature

Qualification Executor

Date



Performance Qualification

Injector Precision

• *Instruments and Fluidics*

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GS50	Dionex	00071255	1.00
LC Module or LC Oven	CC	Dionex	99125564	N/A
Absorbance Detector	AD25	Dionex	00012233	1.02
Autosampler	AS50	Dionex	00052145	1.07

• *Accessories*

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
UV Sample 2	Caffeine, 40 ppm
Eluent A	Water

• *Additional Information*

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
Qualification
Executor/Company: **Mr. Smith / Dionex** Next Qualification: **Feb-01**

• *Test Results Summary*

<i>Test</i>	<i>Result</i>
Injector Precision (Area %RSD)	PASS

Customer Signature

Qualification Executor

Date

• Data for Injector Precision Test

25.0 µl

	Injector Precision
Sample Name	Peak Area AU Caffeine UV_VIS_1
Inj Precision_1	0.049
Inj Precision_2	0.049
Inj Precision_3	0.050
Inj Precision_4	0.050
Inj Precision_5	0.050
Inj Precision_6	0.050
Average:	0.050
Std Dev:	0.000
% RSD:	0.358 %
Limit:	1.000 %
Result:	PASS

Customer Signature

Qualification Executor

Date



Performance Qualification

Detector Linearity

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GS50	Dionex	00071255	1.00
LC Module or LC Oven	CC	Dionex	99125564	N/A
Absorbance Detector	AD25	Dionex	00012233	1.02
Autosampler	AS50	Dionex	00052145	1.07

• Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
UV Sample 1	Caffeine, 15 ppm
UV Sample 2	Caffeine, 40 ppm
UV Sample 3	Caffeine, 60 ppm
UV Sample 4	Caffeine, 80 ppm
UV Sample 5	Caffeine, 100 ppm
Eluent A	Water

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company**

Date: **08-28-00**

Qualification
Executor/Company: **Mr. Smith / Dionex**

Next Qualification Due: **Feb-01**

• Test Results Summary

<i>Test</i>	<i>Result</i>
Detector Linearity (Corr.)	PASS
Detector Linearity (RSD)	PASS

Customer Signature


Qualification Executor

Date

• **Data for Detector Linearity**

Sample Name	Amount ppm Caffeine	Height AU Caffeine UV_VIS_1
Detector linearity_1	15	0.206
Detector linearity_2	40	0.567
Detector linearity_3	60	0.812
Detector linearity_4	80	1.125
Detector linearity_5	100	1.380

• **Linearity Plot**



Cal.Type	Number of Points	Offset	Slope
XLin	5	0.000	0.014

	Correlation Coefficient	% RSD
UV Linearity	1.000	2.208 %
Limit:	0.999	5.000 %
Result:	PASS	PASS

Customer Signature

Qualification Executor

Date



Performance Qualification

Injector Linearity

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GS50	Dionex	00071255	1.00
LC Module or LC Oven	CC	Dionex	99125564	N/A
Absorbance Detector	AD25	Dionex	00012233	1.02
Autosampler	AS50	Dionex	00052145	1.07

• Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
UV Sample 3	Caffeine, 60 ppm
Eluent A	Water

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
Qualification
Executor/Company: **Mr. Smith / Dionex** Next Qualification: **Feb-01**

• Test Results Summary


<i>Test</i>	<i>Result</i>
Injector Linearity (Corr.)	PASS
Injector Linearity (%RSD)	PASS

Customer Signature

Qualification Executor

Date

• **Linearity Plot**



• **Data for Injector Linearity Test**

Sample Name	Retention Time min Caffeine UV_VIS_1	Inj. Vol. µl	Peak Area AU Caffeine UV_VIS_1
Injector linearity_1	0.47	10.0	0.029
Injector linearity_2	0.49	20.0	0.058
Injector linearity_3	0.49	30.0	0.086
Injector linearity_4	0.49	40.0	0.113
Injector linearity_5	0.49	50.0	0.140

Cal. Type	Number of Points	Offset	Slope
LOff	5	0.002	0.003

	Correlation Coefficient	% RSD
Injector Linearity	1.000	0.976 %
Limit:	0.999	1.000 %
Result:	PASS	PASS



Performance Qualification

Injector Carry Over

• *Instruments and Fluidics*

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GS50	Dionex	00071255	1.00
LC Module or LC Oven	CC	Dionex	99125564	N/A
Absorbance Detector	AD25	Dionex	00012233	1.02
Autosampler	AS50	Dionex	00052145	1.07

• *Accessories*

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
UV Sample 6	Caffeine, 1000 ppm
Blank	Water
Eluent A	Water

• *Additional Information*

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
Qualification
Executor/Company: **Mr. Smith / Dionex** Next Qualification: **Feb-01**

• *Test Results Summary*


<i>Test</i>	<i>Result</i>
Carry Over (Area %)	PASS

Customer Signature

Qualification Executor

Date

• **Chromatogram for Carry Over Test**



• **Data for Carry Over Test**

Sample Name	Ret. Time <i>min</i> Caffeine <i>UV_VIS_1</i>	Area <i>AU*min</i> Caffeine <i>UV_VIS_1</i>
Carry over	0.53	0.657
Water	0.51	0.000
Carry over:		0.018 %
Limit:		0.100 %
Result:		PASS

Customer Signature

Qualification Executor

Date



Performance Qualification

Sample Prep Precision

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GS50	Dionex	00071255	1.00
LC Module or LC Oven	CC	Dionex	99125564	N/A
Absorbance Detector	AD25	Dionex	00012233	1.02
Autosampler	AS50	Dionex	00052145	1.07

• Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
UV Sample 6	Caffeine, 1000 ppm
Eluent A	Water
Autosampler Reservoir A	Water

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
Qualification
Executor/Company: **Mr. Smith / Dionex** Next Qualification: **Feb-01**

• Test Results Summary

<i>Test</i>	<i>Result</i>
Sample Prep Precision (%RSD)	PASS

Customer Signature

Qualification Executor

Date

• Data for Sample Prep Precision Test

Sample Name	Peak Area AU Caffeine UV_VIS_1
Sample Prep_1	0.001
Sample Prep_2	0.001
Sample Prep_3	0.001
Average:	0.001
Std Dev:	0.000
% RSD:	0.730 %
Limit:	2.000 %
Result:	PASS

Customer Signature

Qualification Executor

Date



Performance Qualification

Gradient Pump Accuracy and Noise

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GS50	Dionex	00071255	1.00
Absorbance Detector	AD25	Dionex	00012233	1.02

• Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.003" ID PEEK
Eluent A	Water
Eluent B	Water/0.3% acetone (v/v)
Eluent C	Water
Eluent D	Water/0.3% acetone (v/v)

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company**
Qualification
Executor/Company: **Mr. Smith / Dionex**

Date: **08-28-00**
Next Qualification: **Feb-01**

• Test Results Summary


<i>Test</i>	<i>Result</i>
Step Accuracy - A, B	PASS
Step Accuracy - C, D	PASS
Step Noise - A, B	PASS
Step Noise - C, D	PASS

Customer Signature

Qualification Executor

Date

• **Chromatogram of Gradient Pump Accuracy**



• **Step Accuracy - Valves A, B**

Theoretical Value (%)	Absorbance (AU)	Calculated Value (%)	Deviation (%)	Limit (+/- , %)	Result
0.00	0.0000	0.00	0.00	N/A	----
20.00	0.1116	19.67	-0.33	2.00	PASS
50.00	0.2835	49.99	-0.01	2.00	PASS
80.00	0.4520	79.71	-0.29	2.00	PASS
100.00	0.5671	100.00	0.00	2.00	PASS

• **Step Accuracy - Valves C, D**

Theoretical Value (%)	Absorbance (AU)	Calculated Value (%)	Deviation (%)	Limit (+/- , %)	Result
0.00	0.0002	0.00	0.00	N/A	----
20.00	0.1137	20.09	0.09	2.00	PASS
50.00	0.2846	50.37	0.37	2.00	PASS
80.00	0.4524	80.07	0.07	2.00	PASS
100.00	0.5649	100.00	0.00	2.00	PASS

Customer Signature

Qualification Executor

Date

• **Step Noise - Valves A, B**

Step (%)	Noise (AU)	Calculated Noise (%FS)	Limit (%)	Result
0.00	0.0000	0.00%	0.20%	PASS
20.00	0.0002	0.03%	0.20%	PASS
50.00	0.0004	0.06%	0.20%	PASS
80.00	0.0003	0.05%	0.20%	PASS
100.00	0.0002	0.03%	0.20%	PASS

• **Step Noise - Valves C, D**

Step (%)	Noise (AU)	Calculated Noise (%FS)	Limit (%)	Result
0.00	0.0000	0.00%	0.20%	PASS
20.00	0.0002	0.04%	0.20%	PASS
50.00	0.0003	0.05%	0.20%	PASS
80.00	0.0002	0.04%	0.20%	PASS
100.00	0.0002	0.03%	0.20%	PASS

 Customer Signature

 Qualification Executor

 Date



Performance Qualification

AD20 Wavelength Accuracy

• *Instruments and Fluidics*

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Absorbance Detector	AD20	Dionex	97072214	3.02

• *Accessories*

<i>Name</i>	<i>Description</i>
AD20 Holmium Oxide Test Cell	

• *Additional Information*

Customer/Company: **Mr. Turner/ XYZ Company**

Date: **08-28-00**

Qualification

Executor/Company: **Mr. Smith / Dionex**

Next Qualification: **Feb-01**

Customer Signature

Qualification Executor

Date

• Data for the Wavelength Accuracy Test

<i>Wavelength (nm)</i>	<i>Measured Absorbance (AU)</i>
354	0.0000
355	0.0140
356	0.0385
357	0.0621
358	0.0947
359	0.1222
360	0.1356
361	0.1427
362	0.0487
363	0.1470
364	0.1368
365	0.1193
366	0.1070
367	0.0893
368	0.0773

<i>Theoretical Wavelength of Max Absorbance (nm)</i>	<i>Actual Wavelength of Maximum Absorbance (nm)</i>	<i>Limit</i>	<i>Result</i>
361	362.0000	+/- 3nm	PASS

 Customer Signature

 Qualification Executor

 Date



Performance Qualification

Instruments

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GP40	Dionex	98112345	3.40
LC Module or LC Oven	CC	Dionex	99025621	N/A
Chromeleon Datasystem	6.20 Build 527 (Demo-Installation)	Dionex	12345678	N/A
ECD Detector	ED40	Dionex	98103254	3.03
Autosampler	AS50	Dionex	99026688	1.04
Eluent Generator	EG40	Dionex	99037452	2.13

Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
Backpressure Tubing	0.003" ID PEEK, 50"
Blank	Water
ECD Sample 1	Nitrate, 5 ppm
ECD Sample 2	Nitrate, 10.0 ppm
ECD Sample 3	Nitrate, 25.0 pm
ECD Sample 4	Nitrate, 50 ppm
ECD Sample 5	Nitrate, 100 ppm
ECD Sample 6	Nitrate, 1000 ppm
Eluent A	Water
Autosampler Reservoir A	Water
ED/CD Validation Test Cell	
EG40 Validation Test Cell	

Customer Signature

Qualification Executor

Date

Customer signature indicates that all information in the following reports has been reviewed and accepted.

Limits

Test	Customized Limits	Dionex Recommended Limits
Conductivity Noise (nS)	2.000	2.000
Conductivity Drift (nS/hr)	20.000	20.000
Conductivity Linearity (Corr.)	0.999	0.999
Conductivity Linearity (%RSD)	5.000	5.000
DC Amperometry Noise (pA)	4.000	4.000
DC Amperometry Drift (pA/hr)	8.000	8.000
Integrated Amperometry Noise (pC)	4.000	4.000
Integrated Amperometry Drift (pC/hr)	5.000	5.000
Injector Precision (Area %RSD)	1.000	1.000
Carry Over (Area%)	0.100	0.100
Injector Linearity (Corr.)	0.999	0.999
Injector Linearity (%RSD)	1.000	1.000
Sample Prep Precision (%RSD)	2.000	2.000
Standard Bore Gradient Step Accuracy %	2.000	2.000
Standard Bore Pump Noise	0.20%	0.20%
Microbore Gradient Step Accuracy %	2.000	2.000
Microbore Pump Noise	0.50%	0.50%
EG40 Step Current (mA)	See Worksheet	See Worksheet

Additional Information

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
 Qualification Executor/Company: **Mr. Smith / Dionex** Next Qualification: **Feb-01**

Customer Signature

Qualification Executor

Date



Performance Qualification

Conductivity Detector Noise and Drift

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GP40	Dionex	98112345	3.40
ECD Detector	ED40	Dionex	98103254	3.03

• Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
Eluent	Water

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
Qualification
Executor/Company: **Mr. Smith / Dionex** Next Qualification Due: **Feb-01**

• Test Results Summary

<i>Test</i>	<i>Result</i>
Conductivity Noise (nS)	PASS
Conductivity Drift (nS/hr)	PASS

Customer Signature

Qualification Executor

Date

• Data for Detector Noise

Segment No.	Noise (nS)
1	0.376
2	0.438
3	0.391
4	0.486
5	0.394
6	0.402
7	0.427
8	0.365
9	0.309
10	0.324
11	0.286
12	0.291
13	0.296
14	0.296
15	0.318
16	0.331
17	0.328
18	0.299
19	0.288
20	0.341
Average (nS)	0.359
Limit (nS)	2.000
Result	PASS

• Data for Detector Drift


20 Minute Drift, nS	Drift, nS/hr	Limit, nS/hr	Result
0.291	0.873	20.000	PASS

 Customer Signature

 Qualification Executor

 Date

• **Chromatogram of Detector Noise and Drift**



Customer Signature

Qualification Executor

Date



Performance Qualification

Injector Precision

• *Instruments and Fluidics*

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GP40	Dionex	98112345	3.40
Autosampler	AS50	Dionex	99026688	1.04
ECD Detector	ED40	Dionex	98103254	3.03
LC Module or LC Oven	CC	Dionex	99025621	N/A

• *Accessories*

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
ECD Sample 5	Nitrate, 100 ppm
Eluent A	Water

• *Additional Information*

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
Qualification
Executor/Company: **Mr. Smith / Dionex** Next Qualification: **Feb-01**

• *Test Results Summary*

<i>Test</i>	<i>Result</i>
Injector Precision (Area %RSD)	PASS

Customer Signature

Qualification Executor

Date

• **Data for Injector Precision Test**

25.0 µl

<i>Injector Precision</i>	
<i>Sample Name</i>	<i>Peak Area µS Nitrate ECD_1</i>
Inj Precision_1	0.514
Inj Precision_2	0.515
Inj Precision_3	0.513
Inj Precision_4	0.511
Inj Precision_5	0.517
Inj Precision_6	0.516
Average:	0.514
Std Dev:	0.002
% RSD:	0.378 %
Limit:	1.000 %
Result:	PASS

 Customer Signature

 Qualification Executor

 Date



Performance Qualification

Conductivity Detector Linearity

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GP40	Dionex	98112345	3.40
LC Module or LC Oven	CC	Dionex	99025621	N/A
ECD Detector	ED40	Dionex	98103254	3.03
Autosampler	AS50	Dionex	99026688	1.04

• Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
ECD Sample 1	Nitrate, 5 ppm
ECD Sample 2	Nitrate, 10.0 ppm
ECD Sample 3	Nitrate, 25.0 pm
ECD Sample 4	Nitrate, 50 ppm
ECD Sample 5	Nitrate, 100 ppm
Eluent A	Water

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company**

Date: **08-28-00**

Qualification
 Executor/Company: **Mr. Smith / Dionex**

Next Qualification Due: **Feb-01**

• Test Results Summary

<i>Test</i>	<i>Result</i>
Detector Linearity (Corr.)	PASS
Detector Linearity (RSD)	PASS

Customer Signature


Qualification Executor

Date

• **Data for Detector Linearity**

Sample Name	Amount ppm Nitrate	Height µS Nitrate ECD_1
Detector linearity_1	5.0	5.409
Detector linearity_2	10.0	10.337
Detector linearity_3	25.0	24.928
Detector linearity_4	50.0	49.216
Detector linearity_5	100.0	96.009

• **Linearity Plot**



Cal. Type ECD_1	Number of Points ECD_1	Offset ECD_1	Slope ECD_1
XLin	5	0.000	0.978

	Correlation Coefficient	% RSD
ECD Linearity	1.000	4.812 %
Limit:	0.999	5.000 %
Result:	PASS	PASS

Customer Signature

Qualification Executor

Date



Performance Qualification

Injector Linearity

• *Instruments and Fluidics*

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GP40	Dionex	98112345	3.40
Autosampler	AS50	Dionex	99026688	1.04
ECD Detector	ED40	Dionex	98103254	3.03
LC Module or LC Oven	CC	Dionex	99025621	N/A

• *Accessories*

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
ECD Sample 5	Nitrate, 100 ppm
Eluent A	Water

• *Additional Information*

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
Qualification
Executor/Company: **Mr. Smith / Dionex** Next Qualification: **Feb-01**

• *Test Results Summary*


<i>Test</i>	<i>Result</i>
Injector Linearity (Corr.)	PASS
Injector Linearity (%RSD)	PASS

Customer Signature

Qualification Executor

Date

• **Linearity Plot**



• **Data for Injector Linearity Test**

Sample Name	Retention Time min Nitrate ECD_1	Inj. Vol. µl	Peak Area µS Nitrate ECD_1
Injector linearity_1	0.48	10.0	1.005
Injector linearity_2	0.47	20.0	2.005
Injector linearity_3	0.48	30.0	3.003
Injector linearity_4	0.49	40.0	4.020
Injector linearity_5	0.49	50.0	5.014

Cal. Type ECD_1	Number of Points ECD_1	Offset ECD_1	Slope ECD_1
LOff	5	0.000	0.100

	Correlation Coefficient	% RSD
Injector Linearity	1.000	0.191 %
Limit:	0.999	1.000 %
Result:	PASS	PASS

Customer Signature

Qualification Executor

Date



Performance Qualification

Injector Carry Over

• *Instruments and Fluidics*

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GP40	Dionex	98112345	3.40
Autosampler	AS50	Dionex	99026688	1.04
ECD Detector	ED40	Dionex	98103254	3.03
LC Module or LC Oven	CC	Dionex	99025621	N/A

• *Accessories*

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
ECD Sample 6	Nitrate, 1000 ppm
Blank	Water
Eluent A	Water

• *Additional Information*


Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
Qualification
Executor/Company: **Mr. Smith / Dionex** Next Qualification: **Feb-01**

• *Test Results Summary*

<i>Test</i>	<i>Result</i>
Carry Over (Area%)	PASS

Customer Signature _____ Qualification Executor _____ Date _____

• **Chromatogram for Carry Over Test**



• **Data for Carry Over Test**

Sample Name	Ret. Time <i>min</i>	Peak Area <i>µS*min</i>
	Nitrate	Nitrate
	ECD_1	ECD_1
Carry over	0.51	202.975
Water	0.51	0.020
Carry over:		0.010 %
Limit:		0.100 %
Result:		PASS

Customer Signature

Qualification Executor

Date



Performance Qualification

Sample Prep Precision

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GP40	Dionex	98112345	3.40
Autosampler	AS50	Dionex	99026688	1.04
ECD Detector	ED40	Dionex	98103254	3.03
LC Module or LC Oven	CC	Dionex	99025621	N/A

• Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.13 mm (0.005") ID PEEK, 13 m (512")
ECD Sample 6	Nitrate, 1000 ppm
Eluent A	Water
Autosampler Reservoir A	Water

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company** Date: **08-28-00**
 Qualification Executor/Company: **Mr. Smith / Dionex** Next Qualification: **Feb-01**

• Test Results Summary

<i>Test</i>	<i>Result</i>
Sample Prep Precision (%RSD)	PASS

Customer Signature

Qualification Executor

Date

• Data for Sample Prep Precision Test

Sample Name	Peak Area μS Nitrate ECD_1
Sample Prep_1	1.017
Sample Prep_2	1.029
Sample Prep_3	1.020
Average:	1.022
Std Dev:	0.006
% RSD:	0.633 %
Limit:	2.000 %
Result:	PASS

Customer Signature

Qualification Executor

Date



Performance Qualification

Gradient Pump Accuracy and Noise

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Pump	GP40	Dionex	98112345	3.40
ECD Detector	ED40	Dionex	98103254	3.03

• Accessories

<i>Name</i>	<i>Description</i>
Backpressure Tubing	0.003" ID PEEK
Eluent A	Water
Eluent B	1 mM KCL
Eluent C	Water
Eluent D	1 mM KCL

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company**
Qualification
Executor/Company: **Mr. Smith / Dionex**

Date: **08-28-00**
Next Qualification: **Feb-01**

• Test Results Summary


<i>Test</i>	<i>Result</i>
Step Accuracy - A, B	PASS
Step Accuracy - C, D	PASS
Step Noise - A, B	PASS
Step Noise - C, D	PASS

Customer Signature

Qualification Executor

Date

• Chromatogram of Gradient Pump Accuracy



• Step Accuracy - Valves A, B

Theoretical Value (%)	Conductivity (µS)	Calculated Value (%)	Deviation (%)	Limit (+/- , %)	Result
0.00	-0.0061	0.00	0.00	N/A	----
20.00	32.0310	19.60	-0.40	2.00	PASS
50.00	81.8666	50.08	0.08	2.00	PASS
80.00	127.9348	78.26	-1.74	2.00	PASS
100.00	163.4794	100.00	0.00	2.00	PASS

• Step Accuracy - Valves C, D

Theoretical Value (%)	Conductivity (µS)	Calculated Value (%)	Deviation (%)	Limit (+/- , %)	Result
0.00	0.0319	0.00	0.00	N/A	----
20.00	30.4887	18.59	-1.41	2.00	PASS
50.00	80.9523	49.39	-0.61	2.00	PASS
80.00	128.2142	78.23	-1.77	2.00	PASS
100.00	163.8793	100.00	0.00	2.00	PASS

Customer Signature _____

Qualification Executor _____

Date _____

• Step Noise - Valves A, B

Step (%)	Noise (μS)	Calculated Noise (%FS)	Limit (%)	Result
0.00	0.0039	0.00%	0.20%	PASS
20.00	0.1485	0.09%	0.20%	PASS
50.00	0.3392	0.10%	0.20%	PASS
80.00	0.0719	0.04%	0.20%	PASS
100.00	0.0279	0.02%	0.20%	PASS

• Step Noise - Valves C, D

Step (%)	Noise (μS)	Calculated Noise (%FS)	Limit (%)	Result
0.00	0.0038	0.00%	0.20%	PASS
20.00	0.0441	0.03%	0.20%	PASS
50.00	0.0724	0.04%	0.20%	PASS
80.00	0.0858	0.05%	0.20%	PASS
100.00	0.0366	0.02%	0.20%	PASS

 Customer Signature

 Qualification Executor

 Date



Performance Qualification

Integrated Amperometry Detector Noise and Drift

• Instruments and Fluidics

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
ECD Detector	ED40	Dionex	98103254	3.03

• Accessories

Name

ED/CD Validation Test Cell

• Additional Information

Customer/Company: **Mr. Turner/ XYZ Company**

Date: **08-28-00**

Qualification Executor/Company: **Mr. Smith / Dionex**

Next Qualification Due: **Feb-01**

• Test Results Summary

<u>Test</u>	<u>Result</u>
Integrated Amperometry Noise (pC)	PASS
Integrated Amperometry Drift (pC/hr)	PASS

Customer Signature

Qualification Executor

Date

• **Data for Detector Noise**

Segment No.	Noise (pC)
1	0.991
2	1.172
3	1.082
4	1.083
5	1.087
6	0.830
7	1.469
8	1.129
9	1.471
10	0.958
11	0.965
12	0.787
13	0.779
14	0.827
15	0.628
16	1.256
17	1.282
18	1.344
19	0.703
20	1.199
Average (pC)	1.095
Limit (pC)	4.000
Result	PASS

• **Data for Detector Drift**


20 Minute Drift, pC	Drift, pC/hr	Limit, pC/hr	Result
0.023	0.069	5.000	PASS

Customer Signature

Qualification Executor

Date

• **Chromatogram of Detector Noise and Drift**



Customer Signature

Qualification Executor

Date



Performance Qualification

EG40 KOH Current Qualification

• ***Instruments and Fluidics***

<i>Instrument name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>	<i>Moduleware Version</i>
Eluent Generator	EG40	Dionex	99037452	2.13

• ***Accessories***

<i>Name</i>	<i>Description</i>
EG40 Validation Test Cell	
Voltmeter	

• ***Additional Information***

Customer/Company: **Mr. Turner/ XYZ Company**
Qualification
Executor/Company: **Mr. Smith / Dionex**

Date: **08-28-00**
Next Qualification: **Feb-01**

Customer Signature

Qualification Executor

Date

• **Data for KOH Current Qualification Test**

Concentration (mM)	Theoretical Current (mA)	Measured Current (mA)	Limit (+/-)	Result
0	0.00	0.00	0.50	PASS
2	3.22	3.24	0.50	PASS
3	4.82	4.88	0.50	PASS
4	6.43	6.45	0.60	PASS
5	8.04	8.04	0.70	PASS
10	16.08	16.10	1.00	PASS
15	24.12	24.40	1.00	PASS
20	32.16	32.20	1.30	PASS
25	40.20	40.31	1.60	PASS
30	48.24	48.28	1.90	PASS
35	56.28	56.28	2.30	PASS
40	64.32	64.33	2.60	PASS
50	80.41	80.44	3.20	PASS
60	96.49	96.58	3.80	PASS
80	128.65	128.64	5.20	PASS
100	160.81	160.88	6.40	PASS

Customer Signature

Qualification Executor

Date

