



## The TP-Series – the economical line for the Lab



## The new TP-series – simply easy, simply good, simply economic

A new range of compact balances designed to lead the way in providing precision weighing at a great price. Extensive application programs, well-established weighing technology and easy handling makes these balances an ideal solution for standard weighing tasks in the lab, field, industrial or educational environments. Quick accurate weighing powered by state-of-the-art microprocessor technology ensures balance stability and reproducible results.

- Large easy to read high contrast LCD
- Sealed easy to clean keypad
- Stainless Steel weighing pan
- Compact and sturdy design
- Level Indicator and adjustable feet
- Internal, manual calibration in the TP-303i
- Automatic calibration with external weight
- Applications include unit conversion, parts counting, % weighing, display hold, totalizing
- Capable of displaying 14 multiple modes
- Security eyelet
- Bi-directional RS-232 interface
- Weighing Technology "Made in Germany"
- 3 years warranty
- Accessories: Printer, Cal-weights, Dust cover

Model	Capacity	Readability	Pan size	Stability time
TP-214	210 g	0.1 mg	Ø 80 mm	3 sec
TP-303i	310 g	0.001 g	Ø 97 mm	3 sec
TP-303	310 g	0.001 g	Ø 97 mm	3 sec
TP-3002	3100 g	0.01 g	174 x 143 mm	2.5 sec
TP-1502	1500 g	0.01 g	174 x 143 mm	2.5 sec
TP-6101	6100 g	0.1 g	174 x 143 mm	2 sec
TP-3101	3100 g	0.1 g	174 x 143 mm	2 sec
TP-12	12 kg	1 g	174 x 143 mm	1.5 sec



TP-214



TP-303



TP-303i



TP-3002, TP-1502, TP-6101,  
TP-3101, TP-12

### UK office:

Denver House, Sovereign Way  
Trafalgar Business Park  
Downham Market, Norfolk PE38 9SW  
Phone +44.1366.386-242  
Fax +44.1366.386-204

### Denver Instrument GmbH

Robert-Bosch-Breite 10  
37079 Goettingen, Germany  
Phone +49.551.209 77-30  
Fax +49.551.209 77-39  
[www.denverinstrument.com](http://www.denverinstrument.com)  
[info@denverinstrument.de](mailto:info@denverinstrument.de)