

ACQUITY UPLC TUV Detector

The Waters® ACQUITY™ UPLC[™] Tunable UV (TUV) Detector is a tunable, dual-wavelength UV/Visible (UV/Vis) detector specifically designed and optimized for the entire line of ACQUITY UPLC systems. The ACQUITY UPLC TUV Detector's innovative light-guiding flow cell design and High Brightness Lamp are each equipped with Intelligent Chip Technology, allowing for the electronic recording of pertinent component information flow, which is then viewable within the Empower[™] Software database and available for inclusion in a comprehensive status report. Additionally, the design of the ACQUITY UPLC TUV Detector allows for reduced sensitivity to laboratory environmental fluctuations, such as temperature and humidity changes. The ACQUITY UPLC TUV Detector TUV Detector maintains its low-noise electronics and support for sampling rates up to 80 points/s.

OPERATING SPECIFICATIONS¹

Wavelength range	190 to 700 nm	
Wavelength accuracy	±1 nm (via patented ² Erbium filter)	
Linearity range	Deviation at 2.5 AU ≤5%, propylparaben, at 257 nm	
Baseline noise ³	6 x 10 ⁻⁶ AU, at 230 nm, 10 mm shunt cell, 1.0 second filter time constant	
Drift	≤5.0 x 10 ⁻⁴ AU/hour/°C at 230 nm	
Sampling rate	Up to 80 points/s	
Unattended operation	Leak sensors, full diagnostic data captured through console software	

OPTICAL COMPONENT SPECIFICATIONS

Light source	Prealigned, intelligent technology Deuterium lamp (Warranty for one year or 2000 hrs, whichever comes first)	
Dispersion volume	≤0.8 μL (Analytical flow cell) ≤1.5 μL (High sensitivity flow cell)	
Flow cell design	Light-guiding UPLC flow cell with intelligent chip technology	
Path length	10 mm (analytical flow cell)	
Cell volume	500 nL (analytical flow cell)	
Pressure limit	1000 psi (all flow cells)	
Wetted materials	316 stainless steel, fused silica, Teflon™ AF, PEEK (analytical flow cell)	
Optical bandwidth	5 nm	



PHYSICAL/ENVIRONMENTAL SPECIFICATIONS

Dimensions	Width: 34.3 cm (13.5 inches) Height: 19.3 cm (7.6 inches)	
	Depth: 60.9 cm (24.0 inches)	
Weight	12.0 kg (26.5 pounds)	
Operating temperature range	4 to 40 °C (39 to 104 °F)	
Operating humidity range	20% to 80%, non-condensing	
Compliance	CE Mark CSA C-tick and UL	
Audible noise	≤58 dBa	

ELECTRICAL SPECIFICATIONS		
Power requirements	100 to 240 VAC	
Line frequency	50 to 60 Hz	
Power consumption	185 VA (nominal)	
Inputs	One input (inject start)	
Outputs	Two outputs (one event, one analog)	

DETECTOR ORDERING INFORMATION	PART NUMBER
ACQUITY UPLC TUV Detector	176015028
(analytical flow cell included)	
ACQUITY UPLC I-Class Tunable UV Detector	176015092
(low dispersion analytical flow cell included)	
ACQUITY UPLC M-Class Tunable UV Detector (microscale flow cell included)	176016024



OPTIONAL FLOW CELLS

Analytical – 500 nL volume, 10 mm pathlength	205015016
High sensitivity – 2400 nL volume, 25 mm pathlength	205015018
Stainless steel – 1.5 µL volume, 5 mm pathlength	205000610
Titanium – 1.5 μL volume, 5 mm pathlength	205000611
Low dispersion analytical flow cell – 500 nL, 10 mm pathlength	205015033
Low dispersion microbore flow cell – 250 nL, 10 mm pathlength	205015032
Microscale – 100 nL, 10 mm pathlength	205015031

References

- 1. All specifications derive from instrument performance after a 1-hour warm-up time.
- 2. US Patent Numbers: 6,423,249 and 6,783,705.
- 3. ASTM Standard E1657.



Waters, ACQUITY, UPLC, Empower, and The Science of What's Possible are trademarks of Waters Corporation. Teflon is a registered trademark of E.I. du Pont de Nemours and Company. All other trademarks are the property of their respective owners.

©2014 Waters Corporation. Produced in the U.S.A. February 2018 720001322EN TC-KP-PDF

Waters Corporation

34 Maple Street Milford, MA 01757 U.S.A. T: 1 508 478 2000 F: 1 508 872 1990 www.waters.com